Laboratoires Dujardin-Salleron

Catalogue

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Rejoignez notre communauté
Join our community

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Laboratoires Dujardin-Salleron conceive, manufacture and distribute analysis instruments and measurement devices for oenology and spirits. We are pleased to introduce our new catalogue.

Thanks to the close relationship that we maintain with our customers and partners, we were able to develop a wide range of products for many applications. The pages which follow contain more than a thousand references and a lot of information, including details about the proper use of devices.

More than a catalog, this book benefits from our experience to support you throughout your analyzes and measurements. Designed to make your research easier when choosing your laboratory equipment, discover our selection of specialized products:

- *the Dujardin-Salleron's equipment, combining reliability and precision: stills, ebulliometers, titration devices, mustimeters, alcoholmeters,...;*
- *the necessary articles for your day-to-day analyses : laboratory glassware, sampling tools and containers, reagents and standard solutions;*
- *solutions for your hygiene and disinfection processes;*
- *laboratory equipment and furniture;*
- *services of technical training, calibration, maintenance...*

The « Made in France » by Dujardin-Salleron is a criterion of choice, guaranteeing responsiveness from our services and sustainability for your products.

Committed to a process of continuous improvement, we adapt to your needs to offer you total satisfaction.

We hope you find this support useful and we remain at your disposal for any requests or suggestion.

The entire team of Dujardin-Salleron Laboratories

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PH / TITRATION

pH is a crucial parameter in oenology, as it is a good indicator for acidity. It helps to adjust the harvest date and to know the corrections to make during vinification process.

Principe

A solution is systematically defined by its potential hydrogen (pH). It measures the chemical activity of protons (or H^+ hydrogen ions) in solution. It is often associated to acidity or basicity of a solution. In an aqueous solution at 25°C :

- *A solution with $pH < 7$ is defined as acidic ; the lower its pH, the more acidic the solution.*
- *A solution with $pH = 7$ is defined as neutral (water is at $pH = 7$) ;*
- *A solution with $pH > 7$ is defined as basic ; the higher its pH, the more basic the solution.*

Application

Titrimetry is a method used in analytical chemistry to determine the concentration of a chemical in solution. The titration is the implementation of this determination method, involving a chemical reaction between the chemical to measure and a standard reagent that is added progressively. Reaction must be total and measurable (pH, conductivity, color...).

pH metric titrations occur between acids and bases, with an end point detected by pH measurement or color indicator. All acids and bases can be titrated by standard solutions of strong bases or acids.

- *Acidimetric titration : the user drops a base into an acidic volume. Total acidity in musts and wine, which mainly consists of tartaric acid and malic acid, is determined with this method. A standard solution of sodium hydroxide is dropped into a must or wine sample to bring pH to 7. The end point is determined by a color indicator (bromothymol blue) or by a pH meter (DS Titra, Titra EVO) . The volume of sodium hydroxide used allows to calculate the acidity content of the sample.*
- *Alkalimetric titration : the user drops an acid into a basic sample. This method is not used in oenology.*

Titration by redox potential involves another chemical reaction (more details p.14).

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Complete kit for total acidity titration. It includes a 25 ml burette, a stand, a pipette, a crystallizing dish, a reagent and a standard solution for a direct readout.

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Complete kit in plastic case for total acidity titration. It includes an automatic 15 ml burette, a stand, a pipette, a crystallizing dish, a reagent and a standard solution in flexible flask for a direct readout.

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Kit for determination of total acidity. It includes a volumetric graduated flask, a reagent and a standard solution, and a filtering set for musts. Fast and simple operations.

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Dimensions (HxLxP), poids : 40 x 24 x 35 cm, 4 kg

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DS Titra. pH determination, automatic titration of total acidity directly and of volatile acidity after extraction.

Characteristics :

- Automatic titration with a syringe pump activated by a stepper
- Direct display of the result on LCD screen
- Quick and easy-to-use (analysis in less than 1 min)
- Magnetic stirring
- Automatic temperature compensation
- Pre-registered or «customized» titration program
- Ergonomic and aesthetic design
- Excellent repeatability

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Descriptif complet p.16

TITRA EVO. This device enables pH measurement, as well as total acidity determination and volatile acidity after extraction. It avoids the difficult estimation of the detection of color change at the end point, and offers optimum accuracy and reproducibility of the results.

Supplied with accessories and products for around a hundred analysis.

Characteristics :


- Automatic titration
- Optional sample changer 15 positions
- Instinctive touch screen 7"
- Syringe 20 ml
- Memory 4 Go
- Data transfer on flash drive
- Possibility to combine to a feature of titration of SO_2 .

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Apparatus for the extraction of volatile acidity by direct distillation (Duclaux-Gayon method). Laboratory apparatus, with pyrex glassware, simple assembly/dismantling thanks to Rodavis system. Equipped with a rectifying column to stop lactic acid vapors.

116300 Appareil livré sans chauffage

Apparatus delivered without heating accessories

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The same, delivered with butane/propane gas burner

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The following apparatus are designed for the extraction of volatile acidity in conformity with European regulations (Official Journal L272 of 03/10/90) and the OIV standards such as those described in the Recueil des Méthodes Internationales d'Analyse des Vins et des Moûts (Compendium of International Methods of Analysis of Wines and Musts).

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Apparatus for the extraction of volatile acidity by steam distillation (Mothe type). This apparatus with stand is composed of a boiler, an extraction flask, a rectifying column to stop lactic acid vapors and a cooling coil. Extraction time : around 8 minutes.

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Automatic extractor Vola 2000. This apparatus is equipped with a system of instantaneous production of steam (without boiler under pressure) in highly resistant machine-cut aluminium. The high built-in electric heater (2600 W) enables the production of steam by injection, thanks to a dosage pump, of the water necessary for the generator throughout the entire extraction process.

Quick to start and reliable results make the VOLA 2000 advisable for laboratories dealing with wine analysis, spirits, fermented beverages, fruit juices and numerous other agriculture and food industry products (tomato sauces, vinegars, etc...).

The instrument is delivered with all necessary glassware and products for the extraction (except demineralized water).

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Description of burettes in the kits, see box p.7.

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Complete kit for raw volatile acidity titration. It includes a 25 ml burette, a stand, a pipette, a crystallizer, a reagent and a standard solution for a direct readout.

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Complete kit in plastic case for raw volatile acidity titration. It includes an automatic 15 ml burette, a stand, a pipette, a crystallizer, a reagent and a standard solution in flexible flask for a direct readout.

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Complete kit for raw and corrected volatile acidity, including two 25 ml burettes, a stand, a pipette, a crystallizer, an Erlenmeyer flask, standard solutions and reagents.

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Complete kit for acetic acidity titration in vinegars. It includes a graduated burette in acetimetric degrees, a stand, a pipette, a crystallizer, a reagent and a standard solution for a direct readout.

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Complete kit in plastic case for acetic acidity titration in vinegars. It includes an automatic 15 ml burette, a stand, a pipette, a crystallizer, a reagent and a standard solution in flexible flask for a direct readout.

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Complete kit for chromatography on paper. It includes an amber glass ground joint flask, Whatman chromatographic paper, micropipette, reagents and standard samples.

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SO₂ and total acidity kit. It enables to quickly estimate free SO₂ and total acidity. It contains a volumetric graduated flask, a reagent and a standard solution for each titration. Stable solutions, fast and simple operations.

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Complete kit for free and total SO₂ titration (Ripper method). It includes a 25 ml burette, a stand, pipettes, an Erlenmeyer flask, reagents and a standard solution.

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Equipped with an automatic burette, engraved and graduated in mg/l of SO₂ on amber glass flask, with rubber bulb and siphon, for a direct readout.

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Complete kit in plastic case for the titration of free and total SO₂. Free and total SO₂ are successively titrated by a stable standard solution of Iodide/Iodate (KIKIO₃), which progressively gives iodine during titration. All volumes required can be measured with the specific measuring cylinder.

The case contains an automatic 25 ml burette, a stand, a measuring cylinder, an Erlenmeyer flask, reagents and a standard solution in flexible flask.

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Description of burettes included in the kits, see box p.7.

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SO₂ DETERMINATION

SO₂ or sulfurous anhydride, is widely used in oenology. At any step of winemaking process, in particular at the end of alcoholic and malo-lactic fermentations, during ageing and before bottling, it has to be controlled. Its determination is done by titration, directly or after extraction.

By titration

This titration is based on a redox reaction, during which a transfer of electrons occurs. The chemical that loses an electron is called "reducer" and the one that wins one is called "oxidant". Iodine is used as an oxidizing solution to determine SO₂ in acidic conditions.

After all dissolved SO₂ have reacted, iodine is in excess in the solution. This end point is detected by :

- Colorimetry, or Ripper method : a starch solution indicates the end of titration. Added to the sample before the titration, it changes from colorless to blue in the presence of iodine.
- Potentiometry : a double platinum electrode detects precisely the end of titration (DS Sulfi, DS Iodo and TITRA EVO).

In both cases, the volume of iodine used is multiplied by a coefficient, automatically or by the user, to know the content of SO₂ (cf. chart p.121).

Free SO₂ is determined directly, while total SO₂ is determined after releasing combinations of SO₂ with a concentrated solution of sodium hydroxide.

By aspiration-oxidation

The reference method of OIV/EC recommends the use of a Frantz Paul apparatus, characterized by specific dimensions. SO₂ is distilled in acidic conditions by an extraction with gas flow (air or nitrogen). It is collected into an oxidizing solution of hydrogen peroxide (H₂O₂) in which it is converted into sulfuric acid (H₂SO₄), and then quantified by acidimetry.

Free SO₂ is extracted without heating, while total SO₂ is extracted after applying heat to the sample, to enable the extraction of combined SO₂.

This method has the advantage of avoiding the influence of the various compounds of the sample, but it requires a rigorous laboratory practice.

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These automatic titrators determine easily and quickly free and total SO₂. They avoid the difficult estimation of the end of titration with Ripper method. They are particularly adapted to highly colored red wines, lees for distillation, fruit juices, alcoholic drinks and vinegars.

Supplied with the titration glassware and the chemicals needed for approximately 100 analysis.

Characteristics :

- Automatic titration
- Excellent repeatability
- Integrated magnetic stirrer
- Removable stainless steel tray to make the cleaning easier
- Linear guidance of the support of the agitator and the probe
- Power supply : 110-220V - 50/60 Hz - CE

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DS Sulfi. Equipped with auto-zeroing 10 ml burette, graduated 1/20. After the automatic titration, the result is obtained multiplying by a coefficient the volume of solution used (to be read on the burette). Dimensions (HxVxD), weight : 510 x 240 x 350 mm, 2,7 kg

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DS Iodo. Equipped with a syringe pump. After the automatic titration, the result is directly displayed. Dimensions (HxVxD), weight : 400 x 240 x 350 mm, 4 kg

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TITRA EVO. Automatic multiparameter titrator for quick determination of free and total SO₂, pH, volatile acidity on distillate. The device automatically adds the reagents before the titration. It is available with an autosampler to analyse up to 15 samples in a row. Its intuitive touch screen makes it particularly easy to use.

• Iodometric determination of SO₂

With a double platinum electrode that detects an electric current in the sample, the instrument performs the titration and stops automatically when the oxidizing solution is in excess. The device calculates the result depending on the selected programme (free or total SO₂) and displays it directly after a few seconds.

• Determination of pH and acidity

This device allows to determine pH, total acidity directly in the sample, and volatile acidity after extraction, without any need to see a color change at the end of titration. This enables optimum accuracy and reproducibility of the results.

Supplied with titration glassware and the chemicals needed for over a hundred analyses.

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* possible evolution after initial purchase

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Sample changer 15 positions. Supplied separately, it is easy to connect to Titra EVO, to place and to activate, even on a device already set-up.

- Automatic rinsing of the electrode
- Integrated magnetic stirring
- Time per sample : free SO₂ approx.30s, total acidity approx.60s
- Dimensions Titra EVO + changer : 75 x 36 x 44 cm

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2 120720/5 Electrode double de platine (câble non inclus)
Double platinum electrode (without cable)

3 120710/5 Electrode pH (câble non inclus)
pH electrode (without cable)

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Frantz Paul apparatus. Apparatus for the extraction of free and total SO₂ by aspiration-oxydation according to the reference method OIV/CE. The apparatus is composed of a stand, a ground joint extraction flask, a cooler unit, a ground joint bubbling flask and a cylinder. Supplied without heating nor chemicals.

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Electric vacuum mini-compressor. Supplied with hose and adapters for Frantz Paul.

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Heating and fittings for Frantz Paul apparatus. Includes a gas burner, a ring, a gauze and holders.

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Flowmeter for Frantz Paul. 0 to 1 litre/minute. Delivered with connectors.

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Titration kit for Frantz Paul. It enables to directly titrate the sulfuric acid collected in the bubbling flask, with a freshly diluted solution of sodium hydroxide. This kit contains a 10 ml burette, a stand, a standard solution and the accessories for the dilution.

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Fermentest®. Simple and quick method to determine reducing sugars in must and wine. It quickly and easily differentiates dry wines from wines with higher sugar content, by comparison of the color obtained after reaction to a color chart.

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Complete kit for determination of residual sugars in wines and musts with Fehling method. It contains a burette on a stand, a heating system, products for decoloration, Fehling liquor, Erlenmeyer flask and a wooden clamp.

- 196000 Equipé d'une lampe à alcool (chauffe en 6 minutes), d'un entonnoir et de filtres plissés.
Equipped with an alcohol lamp (6-minute heating time), a funnel and folded filters.

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EBULLIOMETRY

Studied and developed by the founders of Laboratoires Dujardin-Salleron, the ebulliometry is a simple and direct method to determine the alcohol content of dry wines and vinegars. It is used worldwide not only by winemakers, but also in other fields of food industry.

Principle

The boiling point is one of the characteristics of a liquid : at a pressure of 1 atmosphere, water boils at 100°C and ethanol at 78,4°C. As wine is a solution mainly composed of water and alcohol, its boiling point depends on its alcohol content. The lower its boiling point, the higher its alcohol content.

Atmospheric pressure and altitude influence the boiling point. Therefore, a calibration is required to determine the environment conditions at the moment of the measurement. So the first step is to measure the boiling point of water, or of a standard wine, to adjust the correspondence table "boiling point / alcohol content". Then the boiling point of the sample can be measured and reported to the correspondence table to know the alcoholic degree.

Sugars in solution have an influence on boiling point. Below 2 g/l, this influence is not important enough to reduce accuracy. Above this threshold, a correction must be applied, or another method must be used.

Recommendations

The measurement of boiling point requires the use of an ebulliometer, traditional or electric. In all cases, a measurement (excluding calibration) lasts about 7-8 minutes. This method is valid for dry wines below 20% A.b.V. with an accuracy of approximately $\pm 0,1\%$ A.b.V.

To obtain more accurate results, some precautions must be taken :

- use the same thermometer or the same probe to perform calibration and measurement ;
- respect scrupulously the volumes of standard solution and sample ;
- work with liquids and an ebulliometer at room temperature ;
- check that refrigeration is efficient to avoid loss of alcohol vapors during heating ;
- make sure there is no sugar in the sample ;
- renew calibration on a regular basis if the ebulliometer is not equipped with a sensor of atmospheric pressure.

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Complete ebulliometer in wooden case, alcohol heating. Developed around 1870, continuously improved and recognised throughout the world, the Dujardin-Salleron Ebulliometer is the simplest instrument for directly and accurately (0,1% AbV) testing the alcoholic content of dry wines with normal composition. Comes complete with wooden case, DS Temp precision thermometer, calculation disk, alcohol lamp, water/wine measurement tube and a box of wicks.

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Special ebulliometer for vinegar testing. The same, with a boiler, the inner surface of which is covered with a protective coating tin/nickel.

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- 3 690114 Tube eau-vin / Water-wine measurement tube
- 4 690112 Boîte de 12 mèches / 12 wicks in a box

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Ebulliometer with DS Temp. With an electronically regulated electric heater and continuous refrigeration, the Ebulliometer enables to accurately determine the alcohol content of dry wines in about 7 minutes. With an ergonomic and original conception, the device is easy to set up and use. It is supplied with all needed accessories to make measurements immediately.

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Ebuliometer. With an electronically regulated electric heater and continuous refrigeration, the Ebuliometer enables to accurately determine the alcohol content of dry wines in about 7 minutes.

Thanks to the automatic atmospheric pressure compensation, there is no need of frequent calibrations.

With an ergonomic and original conception, the device is easy to set up and use. It is supplied with all needed accessories to make measurements immediately.

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Additional methods containing instruments and calculation formula to correct the result of the Ebuliometer according to the density of the sample.

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Ebuliolog. Calculation software to obtain the Alcohol by Volume according to the boiling point and calibration.

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Ebuliotronic. With the same design as the Ebuliometer, this device is equipped with a boiler, the shape and dimensions of which are optimised. Thanks to the integrated calculation tables, the alcohol content is directly displayed in % AbV.

Sample types : dry wines, vinegars, kombucha.

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
Still for plants and liqueurs. Glass apparatus to make tests by steam distillation of plants and infusions, conceived for distillers and perfumers. It works with a heating mantle, not provided (contact us).

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Dujardin-Salleron laboratory still. The Rodaviss system allows a quick and easy dismantling and makes the apparatus airtight. Equipped with a rectifying column and a reflux condenser.

- | | |
|--------|--|
| 140100 | Appareil sans chauffage
<i>Apparatus without heating accessories</i> |
| 140101 | Appareil avec chauffage gaz butane/propane
<i>Apparatus equipped with butane/propane gas burner</i> |

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Mothe laboratory still for all degrees. Equipped with a rectifying bowl, a Graham condenser, and a stopper which ensure a perfect airtightness between the collecting flask and the condenser, it is specially adapted to the analysis of spirits over 25% AbV. Heating parts not included, to be ordered separately.

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Mothe distiller-extractor. It allows a quick steam extraction of alcohol and volatile acids, according to the official method. It is adapted to any kind of wine, liquor (<20% AbV), musts, even if there are sugars or solid pieces. Equipped with a stand and a Graham condenser. Delivered without heating accessories. It is recommended to use a high flow gas burner Ø 40 mm.

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Mothe distiller-extractor with electric heating. Same principle and same applications as the apparatus 140300. Delivered with electric heating 230V - 1000W.

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DISTILLATION - EXTRACTION

Fractional distillation is the type I method recommended by OIV to obtain the Alcohol by Volume (A.b.V.) and volatile acidity, and type IV method for the sorbic acid. When determining the A.b.V., this method is more demanding than ebulliometry, but it is more accurate and it avoids the influence of other compounds on the sample, whatever the sugar and alcohol content.

Principle of laboratory distillation

Distillation is a process used to separate the various liquids in solution : they evaporate one after another, and their vapor is condensed in a condenser to collect a distillate. Laboratory distillation is called "fractional" because of the insertion of a rectifying column in order to achieve a better separation of the liquids of the mixture during extraction.

Applications


Alcohol is extracted by direct heating of sample in traditional glass stills or in automatic devices like Alco 02. This method is recommended for beverages with high alcohol content (from 0 to 80% A.b.V.).

It is also possible to extract by steam the alcohol vapors (recommended method from 0 to 20% A.b.V., possible up to 40% A.b.V. with some conditions). In the same way, traditional models are available along with automatic devices DE EVO and DE 2000. This same method is used for the extraction of volatile acids and sorbic acid.

A third range of instruments includes the traditional and automatic extractors like Vola 2000, exclusively designed for the extraction of volatile acids and sorbic acid by steam.

Determination method after distillation

- Alcohol : density of collected distillate is linked to the A.b.V. It is measured by alcoholmetry, pycnometry, electronic densimeter or hydrostatic balance.
- Volatile acidity : raw volatile acidity is determined by manual or automatic titration. The acidity of free and total SO_2 must be deducted to obtain the corrected volatile acidity.
- Sorbic acid : the concentration of sorbic acid in the distillate is determined by measuring the absorbance at 256 nm by means of a spectrophotometer, comparing it to the ones of a range of standard samples.

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Automatic distillers-extractors designed for the steam distillation of alcohol and the steam extraction of volatile or sorbic acids. All operations comply with European regulation (EEC) Nr.2676/90 and OIV standards described in the Compendium of International Methods of Analysis of Wines and Musts.

Measurement range (alcohol) : 0-40% AbV; from 40 to 80% AbV with 1/2 dilution.

Sample type : wines, spirits, fermented beverages, fruit juices and pastes, tomato sauces, vinegars, etc.

Characteristics :

- Automatic stop combined to the closing of the circulating cooling system
 - Integrated steam generator (2800W)
 - Bubbler with heating rings for an efficient steam extraction
 - Easy drainage of the samples, even if not liquid
 - Equipped with a cleaning gun
 - Easy care and maintenance, easy access to the parts of the device (glass parts, steam generator)
 - Short work time : 5 to 6 min by operation
 - Power supply : 220V - 50 Hz
 - Water supply : minimum pressure 2 bar, recommended temp. 18°C, minimum flow 12 l/min.
 - Dimensions (HxWxD), weight : 93 x 46 x 40 cm - 30 kg
- Supplied with glassware and chemicals needed for the distillation of alcohol and the extraction of volatile acids (except demineralized water).

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DE2000. Selection of the operating time with a timer.

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ALCO 02, automatic direct distiller for spirits and beverages. It complies with the OIV recommendations in the Compendium of Methods of Spirituous Beverages. Equipped with an original pre-heating/heating system, it reduces significantly the time of distillation compared to classical laboratory stills and avoids alcohol losses.

Measuring range : 20-80% AbV.

Sample type : wines, liquors, spirits, etc.

Characteristics :

- Selection of the operating time with a timer
- Automatic stop combined to the closing of the circulating cooling system

- Heating power : 3200 W
 - Constant display of the temperature of the refrigeration liquid
 - Power supply : 220V - 50 Hz
 - Water supply : minimum pressure 2 bar, recommended temp. 18°C, minimum flow 12 l/min.
 - Dimensions (HxWxD), weight : 92 x 44 x 33 cm - 20 kg
- Supplied with glassware and chemicals needed for the distillation of alcohol.*

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DE COMPACT Marcs, automatic distiller by steam extraction. Equipped with high dimensions glass parts and a powerful steam generator, the DE Compact makes the distillation of marc and macerated fruit, in order to determine the alcohol content for upcycling. Complies with French order of August 18, 2014 about declaration and control of the upcycling of winemaking residues. Dimensions (HxWxD), weight : 95 x 35 x 25 cm - 20 kg

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The key points of a laboratory distillation :

- Volume to volume : the volume of collected distillate, after adjustment with distilled water, must be exactly the same as the original sample. The same volumetric flask is used to measure the original sample and collect the distillate. The distillate temperature must be the same ($\pm 2^\circ\text{C}$) as the sample temperature.
- The water running through the condenser must be cold (14-16°C) to ensure an optimal condensation of the alcohol vapors, with a high flow (about 10-12 l/min).

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Mini Free Eau EVO. Designed to optimize the refrigeration of one distiller-extractor. Equipped with a powerful cooling generator, this bath avoids the troubles related to the use of tap water for the refrigeration of automatic distillers (too high temperature, low flow, remote tap, drinkable water waste, cost). It prevents any loss by evaporation, as the first drops of distillate are collected at low temperature.

This compact cool bath is equipped with a thermostatic controller, a water tank and an immersed pump with stirrer. Thanks to its simple structure, it is very easy to install and fill.

Characteristics :

- Nominal cooling efficiency 1700W (Room temperature = 32°C, water at 15°C)
- Dimensions (HxLxP), weight : 55 x 40 x 75 cm, 51 kg (empty)
- Bath temperature : adjustable from -6°C to amb. temperature
- Power supply : 230V-50/60Hz

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
Other baths to cool several distillers simultaneously : Contact us.

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"Chercheur" alcoholmeter 0/100% AbV. Well designed instrument for a perfect vertical stability in the probe and a high accuracy for this kind of instrument. Division : 1% AbV. Length ≈ 300 mm.

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"Chercheur" alcoholmeter calibrated in 2 points, covering the 0/100 % AbV. range with 3 instruments. Division : 0,5 % AbV. Length ≈ 260 mm

147000	0/35
147035	35/70
147070	70/100

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Small size alcoholmeter to measure 100 ml of distillate. Division : 1% AbV. Length ≈ 145 mm.

150000	0/25	151000	0/15
150025	25/40	151010	10/25
150040	40/75		
150075	75/100		

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Precision alcoholmeter calibrated in 3 points, covering the 0/100 % AbV range with 3 instruments. Division : 0,5 % AbV. Length ≈ 300 mm

153000	0/35
153035	35/70
153070	70/100

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Class II precision alcoholmeter calibrated in 3 points, covering the 0/100 % AbV range with 5 instruments. Division : 0,2 % AbV.

Longueur / Length = 310 mm	Longueur / Length = 290 mm
172000 0/20	172060 60/80
172020 20/40	172080 80/100
172040 40/60	172055 55/77

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Class II alcoholmeter calibrated in 3 points. These alcoholmeters comply with the rules defined in French laws. Division : 0,1% AbV. Calibrated at 20°C.

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Standard models, covering the 0/100 % AbV range in 10 instruments.

Longueur / Length = 320 mm	Longueur / Length = 305 mm
174000 0/10	174150 50/60
174010 10/20	174160 60/70
174020 20/30	174170 70/80
174030 30/40	174180 80/90
174040 40/50	174190 90/100

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The same, approved by a certificate of type examination, delivered with a certificate of initial verification to confirm the conformity of the instrument to the reference text, issued after comparison to a COFRAC instrument.

Longueur / Length = 320 mm	Longueur / Length = 305 mm
173000 0/10	173150 50/60
173010 10/20	173160 60/70
173020 20/30	173170 70/80
173030 30/40	173180 80/90
173040 40/50	173190 90/100

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Alcoholmeter kit from 173 series, with an inspected thermometer, in plastic case with protective foam.

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Rack for 6 alcoholmeters. Metal structure with Epoxy coating.

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Bench electronic density meter ALM-155. This densimeter has been especially developed for the enological application, taking advantage of technology and measuring cell of upper-end densimeters, but selecting the sole components for the specific needs of enological field. The Peltier thermostat fixed at 20°C and the specific software enable to offer this densimeter at an incomparable price on the market, with optimal value for money. Sampling by built-in peristaltic pump (approximately 8 ml of sample), simple calibration (air/distilled water), easy cell cleaning, no maintenance needs, auto recording of the results. Measurement range : 0,7 to 1,2 g/cm³, alcohol from 0 to 100,00%. Working temperature : 5 to 35°C. Accuracy : $\pm 5 \times 10^{-5}$ g/cm³, 0,01% AbV.

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Bench electronic density meter DA-850. Measure of the resonant frequency oscillation. Supplied with laptop or tablet. Equipped with Peltier thermostat. The accuracy provided by this instrument is perfectly adapted to oenology and beverages fields. 5 modes of stability according to required accuracy and measurement time, possible expression of density in concentration. Automatic calibration with distilled water or with standard solution. Built-in camera for bubble checking and detection. Minimal sample 1 to 2 ml, injected by a syringe (optional pump). A sample changer can be added. Output : LAN for PC (LIMS connection), USB to download data on flash drive. Output for printer (printer not included). Measurement range : 0 to 3 g/cm³, alcohol from 0 to 100,00%. Working temperature : 0 to 100 °C. Accuracy : $\pm 1 \times 10^{-5}$ g/cm³, 0,01% AbV.

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SPECTROPHOTOMETRY

Spectrophotometry is a quantitative and qualitative analytical method that consists of measuring the absorbance of a solution, i.e. its capacity to absorb light. The absorbance of a sample depends on the concentration of this sample and on the wavelength of the measurement.

Principle

When a light with intensity I_0 goes through a solution, a part of it is absorbed by one or more solute(s). The intensity I of the received light is lower than I_0 . The absorbance is defined by following relation : $A = \log(I_0/I)$. The obtained value does not have any unit. Beer-Lambert law defines the relation between wavelength, the absorbance of a solution, its concentration and the length of the light path.

Application

This principle is directly applied in oenology, to know the color intensity, which is the sum of absorbance values at various wavelengths (420, 520 and 620 nm), and to determine sorbic acid in a distillate obtained by steam extraction, after building a calibration curve with a range of standard samples.


Spectrophotometry is also used to determine easily a wide range of analytical parameters, useful for vinification and ageing, by means of enzymatic and colorimetric kits. The principle is to involve the compound to be determined in an enzymatic or colorimetric reaction that generates a chromophore chemical, which absorbs light at a given wavelength.

Example : to know the malic acid content of a sample, reagents are added to induce the following enzymatic reaction : $L\text{-Malate} + NAD^+ \rightarrow \text{Oxaloacetate} + NADH + H^+$. Then the absorbance at 340 nm is measured, as NADH is chromophore at this wavelength.

We can offer a large range of spectrophotometers in order to meet user's needs :

- analytical laboratories, with references 197200/197205 and 197410/197415, spectrophotometers with touch screen and extended features (spectrum analysis, concentration curves, etc.)
- winemakers who need to control the most common analytical parameters with DS Œnophoto, a compact and light photometer to perform one's own analysis at an affordable cost.


The main analytical parameters that can be measured with enzymatic and colorimetric kits are detailed p.35.

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DS Œenophoto. This photometer with several wave lengths enables fast and accurate analysis by the winemaker for common parameters such as residual sugars, yeast assimilable nitrogen (YAN), malic acid, etc. Completed with usual lab accessories and enzymatic or colorimetric kits, it gives an easy access to these analysis for all wineries.

Characteristics :

- Wavelength : 340 nm and 580 nm
- Absorbance : -0.2 to 3.0 A
- 1 position for macro or semi-micro cuvette with optical path of 10mm
- User-friendly app "DS Œenophoto" to download for free (control and data acquisition)
- Data exportation to an Excel file (csv file)
- Analysis up to 12 samples simultaneously
- Dimensions (HxWxD), weight : 5,5 x 11,5 x 6,5 cm, 230 g

Delivered with 3 micropipettes with disposable tips and a stand, 100 disposables cuvettes and a stand for cuvettes.

To be connected to a smartphone or tablet with Android (with Bluetooth connection), not supplied. Kits to order separately.

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
Spectrophotometers 7200 (Visible) and 7205 (UV/Visible). With diode array technology to produce fast results. Fast and easy set up (color touchscreen). Modes : absorbance, concentration, kinetics, spectrum analysis in 50 points. Small footprint and lightweight. USB ports (communication to a PC). Wide range of accessories available. Equipped with a cuvette position 10x10 mm, and an external cuvette holder.

Spectrophotomètres Spectrophotometers	197200 7200	197205 7205	197410 7410	197415 7415
Gamme spectrale (nm) Wavelength range (nm)	335-800	198-800	320-1000	198-1000
Résolution / Précision Résolution / Accuracy	1 nm / ±2 nm			
Absorbance / Précision Absorbance / Accuracy	-0.300 ... 2.500A ; ± 0.01A			
Bande passante Spectral bandwidth	7 nm	5 nm		
Lampe / Lamp	Tungstène halogène	Xenon	Tungstène	Xenon

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Spectrophotometers 7410 (Visible) and 7415 (UV/Visible). Accurate and reliable resuts. Android interface on a large 7 colour touchscreen. Modes : absorbance, concentration, kinetics, spectrum analysis. Internal memory for data storage (10 Go). USB ports. Multilingual (english, german, spanish, italian). Extensive range of accessories available. Equipped with a single cuvette position 10x10 mm, wide chamber adapted to all accessories.

T.O. 10 mm O.P. 10 mm	Matière (spectre) Material (spectrum)	Taille Size	Quantité Quantity	Usage Use
197812	PS (Vis)	Macro (4,5 ml)	100	Unique Single
197813	PMMA (UV/Vis)	Macro (4,5 ml)	100	Unique Single
197814	PS (Vis)	Semi-micro (2,5 ml)	100	Unique Single
197818	PMMA (UV/Vis)	Semi-micro (2,5 ml)	100	Unique Single
197811	Verre (Vis) Glass (Vis)	Macro (4,5 ml)	1	Lavable Washable
197810	Quartz (UV)	Macro (4,5 ml)	1	Lavable Washable

4 197816 / 197819 Portoir 12 / 24 positions pour cuves.
Cuvette holder with 12 / 24 positions.

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LQA Wine. FTIR analyzer perfectly adapted to routine analysis all along the wine making process, on musts, fermenting musts and finished wines. It combines performance, accuracy, ease-of-use, and short measurement time (less than 45 seconds). Robust and compact instrument with an intuitive user interface, step-by-step workflows, and a built-in connectivity and reporting software, allowing users to master the system with minimal training.

Analytical parameters : see above.

Technical features :

- Equipped with a high performance spectrometer, using fixed paired mirror interferometer technology
- Integrated a measurement cell thermostating system
- Coefficient of variation : $\leq 1\%$
- Analytical capacity : 60 samples per hour
- Sample volume : 5 - 10 ml
- Possible remote connection
- 2 integrated 1L tanks for «zero» and cleaning solutions
- Equipped with a large 12-inch touchscreen, a computer under Windows 10
- Wave Length Range : 400 – 4000 cm⁻¹
- 2 USB ports (for printer, keyboard, mouse, or barcode reader - optional accessories), 1 Ethernet port.

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Thoma: Ruling with 16 big squares of 16 small squares each. Depth 0,1 mm. A big square = 4/1000 mm³.
Malassez: Ruling with 25 big rectangular fields with 20 small squares each. Depth 0,2 mm. A big rectangular field = 1/100 mm³.

- 100136 Simple grille Thoma / Thoma single ruling
- 100137 Double grille Thoma / Thoma double ruling
- 100146 Simple grille Malassez / Malassez single ruling
- 100147 Double grille Malassez / Malassez double ruling

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Cover glass. For counting cell. Box of 10 pieces.

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Slide for microscope, grounded corner 90° (x50)

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Growth media. Specifically conceived for analysis of wines and other fermented drinks. Pre-packaged Petri dishes, available in 55 mm for samples on black membranes with grid, or 90 mm. Aerobic culture at 25°C in incubator (except Lactic Bacteria : anaerobic culture).

Ø 55 mm	Ø 90 mm	Par 10 boîtes de Petri / 10 Petri dishes.
100151	100161	Levures / Yeasts
100157	100167	Levures-moisissures / Yeasts-mold
100152	100162	Brettanomyces / Brettanomyces
100153	100163	Bactéries acétiques / Acetic bacteria
100154	100164	Bactéries lactiques / Lactic bacteria
100155	100165	Bactéries totales / Total bacteria
100156	100166	Flore totale / Total flora

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Anaerobic system. For the incubation of one Petri dish in a bag.

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REFRACTOMETRY

Principe

The refractive index is a characteristic of all liquid - or even pasty - media. The principle is : a light beam that goes from a medium to another changes direction, what determines a refractive angle, specific of the crossed medium. A homogeneous mixture of two components refracts light differently depending on the composition and concentration of this mixture.

The instrument used to determine the refractive angle, the refractometer, is composed of a prism where the sample is dropped and an optical or digital system for direct reading of the result. As the refractive index varies according to the sample temperature, a temperature correction, automatic in most cases, brings the readout to a value at the reference temperature (20°C).

Application

The refractometer is commonly used to determine the sugar content in simple media, such as fruit juices, musts, jam, etc. In this kind of sample, the other substances of dry matter have a very low influence on the refractive angle. The refractive index is converted so that the result is easier to read, generally to Brix degree, or to potential alcohol in oenology.


The refractive index is also used in many fields to control the quality of oils or the concentration of antifreeze liquids for example.

Recommendations

The use of a refractometer is very simple and user-friendly, as the instrument is designed to spread the sample evenly on the prism.

Nevertheless, it is important to take care of a refractometer :

- use preferably a plastic pipette to drop the sample on the prism;
- ensure the device and the sample are at room temperature (difference between both below 5°C);
- check calibration on a frequent basis, by measuring distilled water or a standard solution ; adjust the 0 reading if possible;
- clean the prism surface with water after use and dry with soft tissue;
- avoid shocks or strong vibrations and store away from light and heat.

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Hand held refractometer Master TAL. This ergonomic refractometer is designed to easily and quickly measure the potential alcohol in the must with simplicity and reliability. Besides the automatic distribution of the sample, the instrument is equipped with a high quality of prism, an adjustable eyecup and an automatic temperature compensation from 10 to 30°C. Measuring range : 0/25% AbV PA ±0,2, 0/40% Brix ±0,2.

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RCT 0-22 refractometer. Water resistant and tough model to measure sugar content in grape, it is one of the essential item on the field. Equipped with an automatic temperature compensation, it gives reliable and accurate measures. Supplied with protective case. Range : 0/22% AbV Potential alcohol ±0,2. IP65 protection rating.

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Hand-held refractometer RCT0-22 PA / 0-19 Baumé. With the same design as 196730, it as a double-scale PA/Baumé dedicated to grape musts. Range : 0/22% AbV PA ±0,2, 0-19°Baumé ±0,2.

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Hand-held refractometer Master OE. Same design as Master TAL. It has a double-scale Oechsle/Brix to measure potential alcohol in must. Measuring range: 0/33% Brix - ±0,2%, 0/130° Oechsle - ±1°.

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Digital hand held refractometer PAL-Wine. Compact and water resistant, PAL-Wine is designed to measure the potential alcohol content in grape must. Calibration in accordance with the official method (E.C of 17/09/90). It is equipped with an automatic temperature compensation from 10 to 40°C, with a stainless steel sample stage and a wide liquid crystal display. Measuring range: 0/26% AbV PA ±0,2, 0/40% Brix ±0,2.

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Digital hand held refractometer PAL 84S, 86S, 87S. Same design as PAL-Wine, with specific scales.

- 196713 0/21°Baumé ; $\pm 0.2^\circ$
- 196714 0/240° OeG - 0/53% Brix ; $\pm 1^\circ - \pm 0.2\%$
- 196714B 0/240° Oe - 0/53% Brix ; $\pm 1^\circ - \pm 0.2\%$

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Digital refractometer WM7. It features 7 scales (AP, Brix, Oeschlé, KMW and Baumé.) and is developed to measure quickly and accurately sugar concentration in grape must. It is equipped with an automatic temperature compensation from 5 to 40°C, with a stainless steel sample stage and a wide liquid crystal display.

Gamme de mesure / measuring range :
Alcool probable/potential alcohol 0/26% vol. ± 0.2
Brix 0/45% Brix ± 0.12
Oechsle 0/240° Oe ± 1
KMW (Babo) 0/40° KMW (Babo) ± 0.2
Baumé 0/21° Baumé ± 0.2

2.1.2. Pour fruits, jus de fruits, boissons
For fruit, fruit juice and beverages

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Handheld refractometer, Master-alpha. Same design as Master-TAL. It enables the determination of sugar content in fruits, fruit juice and beverages. Measuring range : 0/32% Brix, $\pm 0.2\%$

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RCT refractometers. Water resistant and tough models to measure sugar content in various sample types. Equipped with an automatic temperature compensation, they give reliable and accurate measures. Supplied with protective case.

	Désignation Designation	Gamme Range	Graduation Graduation
196729	RCT 0-18	0/18% Brix	0,2% Brix
196732	RCT 0-30	0/30% Brix	0,5% Brix
196734	RCT 5-54	5/54% Brix	0,5% Brix

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Digital pocket refractometers. Compact and water resistant, they measure the concentration of many different samples. Automatic temperature compensation from 10 to 75°C.

	Désignation Designation	Gamme Range	Graduation Graduation
196711	PAL-1	0/53% Brix	0,2% Brix
196712	PAL-alpha	0/85% Brix	0,2% Brix

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RX5000i-Plus refractometer. This automatic digital refractometer is equipped with a sample cell on which the sample temperature is adjusted by Peltier effect. It displays quickly and precisely the refractive index, the Brix grade, or the concentration of many liquids. 100 user scales can be defined according to the samples to be analysed. Automatic recording of the 500 latest measurements.
Temperature control range : 5 to 75°C
Measuring range : 1,32420 to 1,58000 nD ± 0.00002 , 0,000 to 100,000 % Brix $\pm 0.001\%$

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pH METRY

pH is an essential tool for winemakers, as it is an indicator of acidity in musts and wine, and it helps to guide vinification. We offer a wide range of pH meters under several categories :

- *pH testers, which give an indication of pH but with low accuracy*
- *portable pH meters to get accurate results on the field*
- *bench pH meters, ideal for use at the laboratory, more ergonomic than portable version to be used on a table with stand for pH electrode*
- *pH meters and conductivity meters to measure conductivity and pH (conductivity meters are also available without pH parameter)*

pH is also measured during titration (more details p.6).

Recommendations for use

When switching on the pH meter, calibrate the electrode with at least two buffer solutions (pH4 and pH7 for example). Buffer solutions must be thrown away approximately one week after the first use. They have a short life after being in contact with air as pH changes quickly and calibration is no longer right. After each measurement, rinse the electrode with distilled water and dry it with a soft paper tissue before performing another analysis or before storage.

With a glass electrode, during a measurement, check that the ceramic junction is properly immersed (1), that the level of liquid in the stem is high enough (2), and that the stopper on the electrode is removed (3). Make sure this stopper is put back on the electrode after use.

Recommendations for storage

To avoid deposits and ensure a fast response, the bulb and the bulb junction with the stem must remain damp. If the pH meter is stored and after each use, clean the electrode and keep it in a storage solution (potassium chloride). Never store the electrode in distilled water. Store the electrode vertically, the bulb at the bottom. A frequent cleaning with a cleaning solution is recommended.

Used at room temperature and in good conditions, an electrode can be used for 1 to 3 years. An electrode not maintained properly or too old increases the response time, gives inaccurate results and requires more and more frequent calibrations.

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Checker electronic pH tester. Low cost electronic digital pH tester. Replaces pH litmus paper. Automatic calibration in 2 points. Range : pH 0 to 14. Accuracy (20°C) : ±0,2.

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GONDO pH tester. This waterproof pH tester is equipped with a digital display showing simultaneously pH and temperature. Automatic calibration in 2 or 3 points. Delivered in case with buffers pH4 and pH7, storage solution, and 4 1,5V batteries. Accuracy : ±0,05.

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PH7 Vio portable pH meter. Professional device with color display, with illumination and contrast. Only 4 keys allow the user to control all functions in an intuitive and easy way. Automatic calibration in 2 or 3 points.

This complete kit includes a pH electrode with integrated temperature sensor, colored pH buffer solutions, carrying case and accessories. Measurement range : pH 0 to 14, resolution pH : 0,1/0,01.

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Five Go portable pH meter. Ergonomic and intuitive, the FiveGo is a user-friendly quality instrument complying with usual requirements of accuracy. It is equipped with an automatic calibration in 3 points, an automatic temperature compensation and data memory for 200 readouts. Delivered with all accessories, ready to use, with a probe holder, a strap, a pH probe with temperature sensor and 2 buffer solution sachets (4.01 and 7.01).

Measurement range : pH 0 to 14. Accuracy : ± 0,01.

INSTRUMENTS DE MESURE

MEASURING INSTRUMENTS

2.2. pH et conductivité (suite) / pH and conductivity (continued)

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Bench pH meter pH50 Vio. This bench pH meter is both easy-to-use and accurate, with a very good value for money. It is equipped with an extra wide LCD with backlight showing simultaneously pH, temperature and measurement stability, and with an automatic temperature compensation. Delivered with combined pH electrode, temperature probe, electrode holder, buffer solutions sachets and electrolyte solution. Measurement range : 0,00 to 14,00 ; Accuracy: $\pm 0,01$.

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Bench pH meter Five Easy. Small foot print, equipped with a practical electrode holder, user-friendly: the FiveEasy is a precise and reliable device, with good value for money. Same technical features and electrode as FiveGo.

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COND7 Vio portable conductivity meter. Professional device with color display, with illumination and contrast. Only 4 keys allow the user to control all functions in an intuitive and easy way. Automatic calibration in 1 to 4 points.

This complete kit includes a probe with integrated temperature sensor, calibration solutions, carrying case and accessories. Measurement range : 0,1 μS to 200 mS ; Accuracy : $\pm 2\%$ f.s.

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Five Go portable conductivity meter. Ergonomic and intuitive, it is equipped with an automatic calibration in 1 point, an automatic temperature compensation and data memory for 200 readouts. Delivered with all accessories, ready to use, with a probe holder, a strap, a probe with temperature sensor and 2 standard solution sachets.

Measurement range : 0,00 $\mu\text{S}/\text{cm}$ to 200 mS/cm ; Accuracy : $\pm 0,5\%$ of displayed value.

INSTRUMENTS DE MESURE

MEASURING INSTRUMENTS

2.2. pH et conductivité (suite) / pH and conductivity (continued)

Mallette du PC70 Vio
PC70 Vio carrying case

Mallette du Seven Go Duo
Seven Go Duo carrying case

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PC5 pH and conductivity tester. Compact and waterproof, it is designed to measure pH, conductivity (EC) and Total Dissolved Solids (TDS) with automatic temperature compensation. Dual level display with 3-color backlight, showing measured value and temperature.

Measurement range : Cond 0-20 mS/cm, TDS : 0-10 ppt, pH -2 to 16, 0-50°C. Accuracy at 20°C : Cond $\pm 1\%$ f.s., pH $\pm 0,1$, $\pm 0,2^\circ\text{C}$.

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PC70 Vio portable pH meter-conductivity meter. Professional device with color display, with illumination and contrast. 6 keys to

control all functions in an intuitive way. USB output (data transfer to PC or power supply). Automatic calibration pH in 1 or 3 points, conductivity 1 to 4 points. Data memory for 500 results. Complete kit in carrying case: 2 electrodes (pH and conductivity with integrated temperature sensor), colored pH buffer solutions, software for PC.

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Seven Go Duo portable pH-conductivity meter. User-friendly 2-way measuring instrument with a high reproductibility level thanks to automatic endpoint recognition, automatic calibration (pH in 3 points) and temperature compensation. Data memory : 99 readouts. Delivered with 2 electrodes with 2-meter cables, rubber case, calibration sachets, sampling flasks, clip for 2 electrodes and a compact carrying case.

Measurement range : 0,10 $\mu\text{S}/\text{cm}$ -500mS/cm, pH 0-14, -5 to 105°C. Accuracy at 20°C : EC $\pm 0,5\%$, pH $\pm 0,01$, $\pm 0,1^\circ\text{C}$

47

INSTRUMENTS DE MESURE

MEASURING INSTRUMENTS

2.2. pH et conductivité (suite) / *pH and conductivity (continued)*

48

2.2.5. Electrodes et sondes / *Probes and electrodes*

Réf. / Ref.	pH	Temp.	Cond.	Prise Connector	Rempliss. Filling	Instrument Instrument
799105	✓					701100
799113	✓	✓				705200
799103	✓	✓		BNC/Cinch	Gel	708100, 758100, 710300
799102	✓	✓		BNC/Cinch	Gel	708200, 758200
799107		✓	✓	BNC/Cinch		718100, 710300
799112		✓	✓	LTW		718200
799111	✓	✓	✓			710110
799106	✓	✓		BNC/Cinch	Gel	710200
799108		✓	✓	LTW		710200

2.2.6. Solutions et accessoires / *Solutions and accessories*

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Electrode holder with flexible arm.

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Cleaning solution for pH electrode, 250 ml.

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Storage solution for pH electrode, 500 ml.

5 971050 Solution électrolyte 30 ml / Electrolyte solution 30 ml

6 Solutions tampon / *Buffer solutions*

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
100 g range hydrometer, divided in 1 g/dm³. Length ≈ 280 mm

330007	700/800	330012	1200/1300
330008	800/900	330013	1300/1400
330009	900/1000	330014	1400/1500
330095	950/1050	330015	1500/1600
330010	1000/1100	330016	1600/1700
330105	1050/1150	330017	1700/1800
330011	1100/1200		

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Precision hydrometer : these hydrometers, covering 30 grams per instrument, comply with the rules defined in ISO norm. Division : 0,2 g/dm³. Calibrated at 20°C.

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Mustimeter with double scale (density 980-1130 g/dm³ and potential alcohol). According to EC rule Nr.2676/90. It allows to measure the sugar content in grape must according to its density and the potential alcohol content, to be obtained after complete fermentation.

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Guyot glucometer for grape musts, with triple scale (Baumé degree, potential alcohol and sugars in kg/hl). Length ≈ 300 mm

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Cider hydrometer. 1000/1100 hydrometer for apple must and ciders with double scale (density and potential alcohol).

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Beer hydrometer. Double scale (density 1000-1100 g/dm³, Plato degree) to check fermentation progress.

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Use of mustimeter : handling precautions

- The instrument must be perfectly clean and dry : clean with a solution of 20% sodium hydroxide, then rinse with distilled water and dry it with a clean tissue.

- Use it in a glass cylinder without graduation, straight (not conical). Do not use the instrument directly in the tipping trailer or tank.


- Make sure that the sample measured is representative. Remove rough particles to avoid wrong measurements (by settling, sieving, etc.).

- Immerse slowly the instrument in the sample, holding it at the top of the stem, giving an up and down movement, then release it with a rotary movement.

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Melli-mustimeter with triple scale (density, weight of honey in kg/hl and potential alcohol for honey wine).

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Hydrometers for various applications

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«Tirage» hydrometer 990/1005 g/dm³ divided in 0,1. Precise and sensitive model designed to be used for «tirage» wines. Calibrated at 15°C to be used in wine cellars. Length ≈ 325 mm.

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Extracto-enometer. Precision hydrometer to determine easily the dry extract in wines with hydrometric method (983-1003 g/dm³). Delivered with instructions and tables.

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Divided in 0,2 g/dm³, with an inspection certificate to confirm the conformity of the instrument to the reference text, issued after comparison to a COFRAC standard instrument. Length ≈ 325 mm

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Salinometer. Precision hydrometer to determine seawater density.

339005 1000/1030 divisé en 1 g/dm³, long. 200 mm
1000/1030 divided in 1 g/dm³, length 200 mm

339000 1005/1030 divisé en 1/5 g/dm³, long. 300 mm
1005/1030 divided in 1/5 g/dm³, length 300 mm

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Portable electronic densimeter DA130N. This easy handling instrument, based on the measure of the resonance frequency, allows the measurement of density of many liquids, even viscous (wines, spirits, etc.). Automatic temperature compensation from 0 to 40°C. Output : IR port (Wire not included). Weight : 360 g
Measure range : 0 to 2 g/cm³
Accuracy : ± 0,001 g/cm³ - 0,1°C

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Portable electronic densimeter Densito. With excellent ergonomics and an intuitive programme guided by a joystick, Densito offers fluid and fast process, with accuracy and repeatability. Automatic temperature compensation from 0 to 50°C. Output : USB port. Weight : 355 g
Measure range : 0 to 3 g/cm³ - Accuracy : ± 0,001 g/cm³ - 0,2°C

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The same, "Pro" version. Equipped with a RFID and bar-code to improve traceability of measurements.

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Red liquid filled laboratory thermometer

	Portée en °C Range in °C	Division Division
403009	-10/+60	1
403010	-10/+110	1
403007	-5/+35	1/2
403020	0/+35	1/2
403008	-5/+50	1/2
403021	0/+50	1/2
403022	-10/+50	1/2
403011	-10/+60	1/2
403023	0/+35	1/10

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PTFE protective sleeve. Application of a PTFE sleeve on a stem thermometer to protect marking from alcohols and to reduce the risks in case of breakage.

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Inspected thermometer with PTFE protective sleeve.

	Portée en °C Range in °C	Division Division
406020	0/+35	1/2
406021	0/+50	1/2
406023	0/+35	1/10

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Enclosed scale thermometer, red liquid, 0/+40°C, division 1/2.

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Red spirit filled thermometer with protective sleeve

	Portée en °C / Range in °C	Gaine / Sleeve
408000	-10/+60	Polypro
408003	-10/+120	Polypro
409000	-10/+120	Bois / Wood

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Use of an expansion thermometer

Before use, check that the liquid in the column is not separated due to rough handling in transport or incorrect storage. In case of separated liquid, it is necessary to reunite it by heating, cooling or centrifuging. Position the thermometer vertically immediately after this operation. Check that the thermometer is used in its measuring range. Store vertically.

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DSTemp. Digital Thermometer with high accuracy -20/+120°C with Pt1000 probe. Resolution : 0,01°C.

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Electronic thermometer for tank. K type thermocouple thermometer with LCD, -50/+300°C, ±1°C. HOLD function.

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Stainless steel probes for digital thermometer 410103 and 411040.

- 411039 30 cm, rigide / 30 cm, rigid
411101 150 cm, rigide / 150 cm, rigid
411102 210 cm, rigide / 210 cm, rigid
411103 150 cm, souple / 150 cm, flexible

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HQ1130 portable oxymeter. Equipped with a LDO probe : no membranes to clean or replace, no electrolyte solution to refill, reduced calibration frequency, adapted to any kind of wine or beverage. The only maintenance is the replacement of the sensor cap only once a year. Automatic atmospheric pressure compensation. Data memory : 5000 readouts. Delivered with carrying case.

*Measurement range : 1 to 200 % or 0,1 to 20,0 mg/l
Accuracy : ±0,1 from 0 to 8 mg/L, ±0,2 above 8 mg/L*

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Bentotest, complete kit for the quick determination of proteic stability. It includes a 25 ml cylinder, reagents, a funnel and folded filters.

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Portable turbidimeter HI98703-02. Instrument with good value for money to determine turbidity during white wines vinification. Nephelometric method ratio of scattered and transmitted light. 3 measuring modes (normal, average, and continuous). Tungsten lamp. Data memory : 200 records. RS232 and USB outputs. Delivered in carrying case with cells, silicon oil, calibration set, 12V adapter and batteries.

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Portable turbidimeter 2100Q. Reference instrument, user-friendly, to determine turbidity at all steps of vinification and before bottling. Ratio turbidimetric determination using a primary nephelometric light scatter signal (90°) and transmitted light scatter signal. The two-detector optical system compensates color in the sample, light fluctuation, and stray light.

Tungsten lamp. Data memory : 500 records.

Delivered in carrying case with cells, calibration kit STABLCAL, silicon oil. Measurement range : 0 to 1000 NTU (Accuracy $\pm 2\%$)

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Bench turbidimeter TL2300. Covering a wide range of turbidity, it performs accurate measurements based on a ratio of signals from 4 sensors. Large full color touch screen for setup, calibration and measurements easier. Tungsten lamp. USB-A Port. Delivered with sample cells, StablCal and Gelex standard samples, silicone oil, and dust cover. Measurement range : 0 to 4000 NTU Accuracy from 0 to 1000 NTU : $\pm 2\%$ + 0,01 NTU

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Checkstab iProtein. The semi-automatic analyzer, iProtein, is the first complete system for the Heat Test for Stability, equipped with a software for the evaluation of Protein Stability. This innovative system has been designed to simplify, computerize and trace this analytical parameter considered in the oenological field. iProtein is a semi-automatic analyzer, equipped with a thermostatic bath for heating, a cryostat for cooling and an integrated turbidimeter adapted to this type of analysis. Simultaneous analysis of 15 samples, identified by RFID technology. Customizable methods and automatic calculation of needed bentonite for wine stabilization.

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The clogging index (CI) gives information about the wine filterability. Its determination requires a 2-bar pressure on 400 ml of wine filtered on a membrane, Ø25 mm, porosity 0,65 µm. It is based on the measurement of the time of filtration of successive parts of a wine volume at room temperature. A stainless tank with lid on a stand is mounted on the filter holder, and it is completed with a manometer.

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DISSOLVED CO₂

All wines, still and sparkling, contain CO₂ (carbonic gas). In particular, CO₂ plays an important part in still wine tasting. Vinification must take into account this parameter to assure that the CO₂ content during bottling is adjusted to the wine profile.

The determination of dissolved CO₂ in still wines is done easily with a Carbodoseur, giving quite accurate results if the measurement is performed in appropriate conditions, in particular at low temperature (sample below 10°C). Some more sophisticated devices enable to obtain fine-tuned values thanks to NDIR (Non Dispersive Infra Red) technology.

PRESSURE

In bottles of sparkling and semi-sparkling wines, excess pressure due to carbon dioxide can be measured with an aphrometer (pressure gauge). Their structure may vary depending on the closure type of the bottle, but most of them include a pressure gauge (manometer), a holding system and a piercing system.

The measurement must be done in the following way: the aphrometer is placed on the bottle and the fixing bracket positioned so that the instrument is stable on the bottle and there is airtightness between the inner space of the bottle and the aphrometer body. After piercing the crown cap or the cork with the aphrometer, the bottle must be strongly shaken until pressure is stable. Then the reading is done. In the case of permanent aphrometer, the device replaces the crown cap and is placed directly on the bottle ring. It allows to control the pressure increase all along the secondary fermentation.

The measurement must be done on bottles whose temperature has been stable for at least 24 hours. If the measurement temperature is different from 20°C, a correction must be applied, multiplying the measured pressure by the appropriate coefficient.

The aphrometer can also be used on bottles filled with still wines and closed with a classic cork, in order to check bottling quality. It is available for this use under a simple structure, without any holding system, and with a vacuum/pressure gauge. Bottles with screw cap require an aphrometer with fixing bracket.

To comply with OIV norm "Overpressure measurement of sparkling wines", the aphrometers must be regularly verified (at least on an annual basis) by means of a calibration bench.

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Carbodoseur : this is the simplest and most appropriate method to control dissolved CO₂ at the winehouse. The measure is done with a stoppered cylinder filled with the sample to analyze. Once shaken, the cylinder spills out a certain volume of the sample, proportionately with the CO₂ content. Average accuracy of 50 mg/l. Comes complete in a plastic case with calculation tables.

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Simple aphrometers to control the pressure on bottles stoppered with cork at the level of the neck, for still or slightly sparkling wines.

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INSTRUMENTS DE MESURE

MEASURING INSTRUMENTS

2.1. Pression (suite) / Pressure (continued)

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Permanent aphrometer with manometer 0/10 bar. Dedicated to permanent control of pressure increase during fermentation in bottle, it must be placed on the neck (Ø29 mm) directly after bottling, without any cap or cork. Accuracy is guaranteed by a perfect airtightness of holding system.

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"Cap" aphrometers. Equipped with a 0/10 bar manometer, a piercing system and a fixing system. The fixing ring is removable and replaceable. A gasket ensures airtightness before piercing the crown cap.

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«Cork» aphrometer for sparkling wines. Equipped with a 0/10 bar manometer, it allows to measure pressure instantaneously on bottles stoppered with cork. The detachable tip falls down inside the bottle after going through the cork to allow pressure measurement by the end of the needle.

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Options (integrated bidule, synthetic cork): contact us.

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CheckStab. Quick and easy determination of tartaric stability based on electrochemical conductivity measurements under various conditions.

Different methods are used for tartaric stability measurement :

- Mini contact method determines KHT precipitation, expressed in µSiemens. This value indicates the stability of the wine being tested.
- Saturation point in °C, automatically registered for each tested sample, without any intervention by the operator.
- Mini contact with Forecast method determines KHT precipitation expressed in µSiemens. By using CheckFuture®, it predicts precipitation of crystals.

The Check Stab also allows to conduct calcium saturation point analysis. In this case, Calcium L-Tartrate Hydrate C4H4CaO6 is needed (supplied as an option).

Advantages :

- Check control function allows easy and fast analyses.
- 7,1" color touch screen display
- Silent ventilation system
- Automatic KHT dispensing
- USB interface to download data and upgrade the programme
- Standby function for energy saving

*Supplied with analytical software system to manage data
 Dimensions (HxLxD), weight : 505 x 280 x 420 mm, 28 kg*

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Spark M20. Thanks to an innovative air mixing process in the sample, it allows distillers and control services to easily and quickly verify the quality of spirits. Indeed, these products can be contaminated with methanol, which can cause blindness and lead to death. Portable device with intuitive interface. Analysis in 2 minutes. Measurement range : 0-10 % vol., resolution 0,01%, accuracy 10% of value.

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Spark M30. The same, with settings dedicated to detection of methanol traces in products with low alcohol content, after addition of DMDC.

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Mini Logger USB is a single-use temperature recorder, in a waterproof pouch, with s/n, programmed to measure temperature on a defined period, to control the temperature conditions during storage or shipment. PDF report. LEDs on the key give information about overrunning the low and high temperature limits. Measurement range : -40 to +80°C; Accuracy : $\pm 0,5^\circ\text{C}$.

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Hair tension hydrometer. Hygrometric scale 0 - 100%. Ø100 mm.

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Digital hygrometer-thermometer. Simultaneous display of hygrometry and temperature. MAX/MIN function.

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Data logger for temperature and humidity KH50. Perfectly adapted to winehouses, this data logger is very simple to set up and read. Large LCD. 2 configurable setpoint alarms, with visual signal (LED). Magnetic mounting. Up to 16000 recording points. Measurement range : 0 - 100 % RH / -20 to +70 °C

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LABORATORY GLASSWARE AND KITS

3.1. Burettes et accessoires / Burettes and accessories

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Head for 24/29 ground joint flask for burette and bulb

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Schilling burette for a quick titration, automatic, mounted on a flexible PE flask and stand, equipped with pinch cock with titration screw for precise regulation of the flow.

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Burette with stopcock, engraved in mg/l of SO₂

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LABORATORY GLASSWARE AND KITS

3.1. Burettes et accessoires (suite) / Burettes and accessories (continued)

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Burette with stopcock, engraved in acetimetric degrees

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- Pinces et supports / Clamps and support

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25 ml Brand bottle-top digital burette, piston-type, easy-to-use : rapid filling of the burette and precise titrations with the same hand wheels, 360° directional system. Easy and quick calibration. Quick air purge without air bubbles. Resolution 0,01 ml - Accuracy <0,1% Delivered without flask, with an injection tube, a telescopic filling tube, an assembly key, PP adapters and 2 AAA batteries.

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2,5 l amber glass flask for BRAND bottle-top burette

3.2. Cristallisoirs / Crystallizing dishes

3 Cristallisoir à bec en verre. DIN 12338.

Crystallizing dish with pouring spout, DIN 12338.

	Diamètre Diameter	Capacité Capacity	Hauteur Height
627040	40 mm	20 ml	25 mm
627050	50 mm	40 ml	30 mm
627060	60 mm	60 ml	35 mm
627070	70 mm	100 ml	40 mm
627080	80 mm	150 ml	45 mm
627090	95 mm	300 ml	55 mm
627100	115 mm	400 ml	65 mm
627150	140 mm	900 ml	75 mm
627170	190 mm	2000 ml	90 mm
627190	230 mm	3500 ml	100 mm

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Crystallizing dish without pouring spout, DIN 12337. Diameter 70 mm, volume 100 ml, height 40 mm.

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Graduated glass beaker, with pouring spout, low form, ISO3819.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
613050	50 ml	42 mm	60 mm
613100	100 ml	50 mm	70 mm
613150	150 ml	60 mm	80 mm
613250	250 ml	70 mm	95 mm
613400	400 ml	80 mm	110 mm
613600	600 ml	90 mm	125 mm
613800	800 ml	100 mm	135 mm
613901	1000 ml	105 mm	145 mm
613902	2000 ml	130 mm	185 mm

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Graduated PP beaker, with pouring spout, low form, blue marking, ISO7056.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
613050P	50 ml	41 mm	60 mm
613100P	100 ml	51 mm	72 mm
613250P	250 ml	71 mm	95 mm
613500P	500 ml	87 mm	119 mm
613901P	1000 ml	109 mm	147 mm

3 Bécher gradué en verre. Bécher sans bec forme haute. Graduated glass beaker, without pouring spout, high form.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
613050H	50 ml	38 mm	70 mm
613150H	150 ml	54 mm	95 mm

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Graduated glass Erlenmeyer flask, narrow neck, ISO1773

	Capacité Capacity	Diamètre base/col Base/neck diameter	Hauteur Height
615100E	100 ml	64/22 mm	105 mm
615250E	250 ml	85/34 mm	140 mm
615500E	500 ml	105/34 mm	175 mm
615901E	1000 ml	131/42 mm	220 mm
615902E	2000 ml	166/50 mm	280 mm

5 Fiole Erlenmeyer col large en verre. Fiole graduée, DIN12385.

Graduated glass Erlenmeyer flask, wide neck, DIN12385

	Capacité Capacity	Diamètre base/col Base/neck diameter	Hauteur Height
615100L	100 ml	64/34 mm	105 mm
615250L	250 ml	85/50 mm	135 mm
615500L	500 ml	105/50 mm	170 mm
615901L	1000 ml	131/50 mm	280 mm

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Non graduated glass cylinders, tall form without pouring spout, especially designed for hydrometers and alcoholmeters.

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450 ml cylinder for mustimeters, with thermometer slot.

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325 ml calibrated cylinder for high degree alcoholmeters (173 and 174 series, above 50% AbV). Depth : 310 mm - int. Ø : 36 mm

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400 ml calibrated cylinder for low degree alcoholmeters (173 and 174 series, below 50% AbV). Depth : 330 mm - int. Ø : 39 mm

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450 ml calibrated cylinder for all alcoholmeters from 173 and 174 series.

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LABORATORY GLASSWARE AND KITS

3.5. Eprouvettes (suite) / Cylinders (continued)

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Graduated glass cylinder, class A, tall form, with pouring spout and hexagonal base, ISO4788.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
623005	5 ml	12,9 mm	115 mm
623010	10 ml	15,5 mm	135 mm
623025	25 ml	21,3 mm	160 mm
623050	50 ml	26,0 mm	195 mm
623100	100 ml	31,3 mm	240 mm
623250	250 ml	40,7 mm	320 mm
623500	500 ml	53,2 mm	380 mm
623901	1000 ml	67,0 mm	465 mm
623902	2000 ml	85,0 mm	535 mm

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Graduated glass cylinder with plastic stopper, class A, tall form, blue graduation and hexagonal base, ISO4788.

	Capacité Capacity	Diamètre Diameter	Hauteur Height	Rodage Ground joint
625100	100 ml	31,3 mm	240 mm	24/29
625250	250 ml	40,7 mm	320 mm	29/32
625500	500 ml	53,2 mm	380 mm	34/35
625901	1000 ml	67,0 mm	465 mm	45/40
625902	2000 ml	85,0 mm	535 mm	45/40

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Graduated translucent polypropylene (PP) cylinder, tall form, with moulded graduations, pentagonal base, and pouring spout, in accordance with ISO6706.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
624050R	50 ml	25,5 mm	199 mm
624100R	100 ml	30,5 mm	249 mm
624250R	250 ml	41,5 mm	315 mm
624500R	500 ml	55,0 mm	361 mm
624901R	1000 ml	66,0 mm	439 mm
624902R	2000 ml	84,0 mm	531 mm

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Graduated PMP (polymethylpentene) cylinder, crystal clear, tall form, with blue graduations, pentagonal base, and pouring spout, in accordance with ISO6706.

	Capacité Capacity	Diamètre Diameter	Hauteur Height
626050	50 ml	26,0 mm	199 mm
626100	100 ml	31,0 mm	249 mm
626250	250 ml	41,5 mm	315 mm
626500	500 ml	55,0 mm	361 mm
626901	1000 ml	66,0 mm	438 mm
626902	2000 ml	84 mm	531 mm

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Volumetric glass pipette, Class AS, with safety bulb, ISO648.

1 trait 1 mark	2 traits 2 marks	Capacité Capacity	Hauteur Height
602001	603001	1 ml	310 mm
602002	603002	2 ml	330 mm
	603004	4 ml	350 mm
602005	603005	5 ml	400 mm
602010	603010	10 ml	440 mm
602015	603015	15 ml	455 mm
602020	603020	20 ml	510 mm
602025	603025	25 ml	520 mm
602050	603050	50 ml	540 mm

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Graduated glass pipette, 1/10, Class AS, with safety bulb, ISO835.

	Capacité Capacity	Hauteur Height
604001	1 ml	350 mm
604002	2 ml	350 mm
604005	5 ml	360 mm
604010	10 ml	360 mm
604020	20 ml	450 mm
604025	25 ml	450 mm
604050	50 ml	720 mm

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Pasteur pipette. Disposable PE pipette with bulb and graduations. Capacity 3 ml, graduation 0,5 ml.

605103	Stérile, en emballage individuel Sterile, individually wrapped
605203	Non stérile / non-sterile

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Calibrex dispensers with variable volume. With high tech materials and ergonomic design, this range is dedicated to precise dispensing of liquids from a large variety of bottles and flasks, with long lasting performance. Instant volume setting. Autoclavable at 121°C fully assembled. PFA glass coated plunger that enables trouble free distribution of salt solutions, weak and strong acids, as well as bases. Supplied without flask, with 3 PP bottle adapters

	Volume Volume	Divisions Division	Précision / Accuracy Vol.min. Vol.max.	
600161	1 - 5 ml	0,1 ml	< 0,5%	< 0,1%
600162	1 - 10 ml	0,25 ml	< 0,5%	< 0,1%
600163	2,5 - 25 ml	0,5 ml	< 0,5%	< 0,1%
600164	5 - 50 ml	1,0 ml	< 0,5%	< 0,1%

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Volumetric glass flask, class A, with standard ground joint and PET stopper, ISO1042, DIN12664. Delivered with batch certificate.

	Capacité Capacity	Diamètre Diameter	Hauteur Height	Rodage RIN Stopper size
618050D	50 ml	50 mm	140 mm	12/21
618100D	100 ml	60 mm	170 mm	14/23
618200D	200 ml	75 mm	210 mm	14/23
618250D	250 ml	80 mm	220 mm	14/23
618500D	500 ml	100 mm	260 mm	19/26
618901D	1000 ml	125 mm	300 mm	24/29
618902D	2000 ml	160 mm	370 mm	29/32

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Volumetric glass flask, with one graduation mark, class A, especially adapted for laboratory stills.

	Capacité / Capacity	Hauteur / Height
619200	200 ml	180 mm
619250	250 ml	185 mm
619300	300 ml	200 mm

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Glass flask with 2 marks for sugar analysis, without stopper.

	Capacité / Capacity	Hauteur / Height
616050	50/55 ml	155 mm
616100	100/110 ml	180 mm
616200	200/220 ml	200 mm

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LABORATORY GLASSWARE AND KITS

3.10. Ballons et accessoires / Flasks and accessories

Rodage femelle et mâle a/b
Female and male ground joint a/b

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Very precise volumetric flask to check fill volumes. Factory calibration certificate not included, available upon order.

	Capacité Capacity	Graduation Graduation	Hauteur Height
617037	375 ml	365-385 ml	340 mm
617075	750 ml	680-780 ml	460 mm
617100	1000 ml	940-1040 ml	490 mm
617150	1500 ml	1440-1540 ml	520 mm

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Wooden protective case for calibrated fill flasks.

	Pour fiole For flask
617037C	617037 (375 ml)
617075C	617075 (750 ml)
617100C	617100 (1000 ml)
617150C	617150 (1500 ml)

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Automatic wooden tong. Length 250 mm. Grips a minimum diameter of 10 mm, maximum 35 mm.

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Round bottom glass flask, ground joint neck (Rodaviss).

	Capacité Capacity	Ø ext. corps Body outer Ø	Rodage Ground joint	Hauteur Height
610250S	250 ml	85 mm	29/32	135 mm
610500S	500 ml	105 mm	29/32	150 mm
610901S	1000 ml	131 mm	29/32	200 mm

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Other flasks with special neck : contact us.

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PP flask holder for round bottom flasks up to 10 litres. Ø160 mm. It supports hot flasks up to 120°C. Autoclavable.

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- Pince bois / Wooden tong
- Chauffe-ballons / Heating mantle

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Glass funnels with short stem, for filter paper, 60° angle, ISO4798.

	Capacité Capacity	Ø supérieur Upper Ø	Ø tige Stem Ø	Long.tige Stem length
363001	8 ml	35 mm	6 mm	35 mm
363002	15 ml	45 mm	6 mm	45 mm
363003	30 ml	55 mm	8 mm	55 mm
363004	115 ml	75 mm	8 mm	75 mm
363005	140 ml	85 mm	9 mm	85 mm
363006	190 ml	100 mm	10 mm	100 mm

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PP funnels with short stem, 60° angle, handle with eyelet.

	Capacité Capacity	Ø supérieur Upper Ø	Ø tige Stem Ø	Long.tige Stem length
363010	88 ml	75 mm	7 mm	55 mm
363011	222 ml	100 mm	8 mm	77 mm
363012	342 ml	120 mm	11 mm	90 mm
363013	817 ml	150 mm	15 mm	95 mm

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Double funnel holder. PP holder with epoxy base and stainless steel rod. For funnels with diameter 50 to 120 mm.

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Standard filter paper, smooth, medium filtration for general use in laboratory. Characteristics : 64 g/m², thickness 0,15 mm, porosity 10-20 µm. Pack of 100 folded filters.

	N°	Diamètre Diameter	Pour entonnoir réf. For funnel ref.
650300	00B	100 mm	363001, 363002
650310	0B	130 mm	363003
650311	1B	150 mm	363004, 363010
650312	2B	190 mm	363005
650313	3B	250 mm	363006, 363011, 363012
650314	4B	330 mm	363013
650315	5B	400 mm	
650316	6B	450 mm	
650317	7B	500 mm	

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Quick filtration system. Complete kit : 10 ml glass syringe with Luerlock connector, a filter holder (Ø25 mm) and 100 filter disks.

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
Filter disks (x100) for quick filtration system. Ø25 mm, por. 10-20 µm.

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LABORATORY GLASSWARE AND KITS

3.12. Matériel de filtration (suite) / Filtration equipment (continued)

Mini-compresseur non inclus
Mini-compressor not included

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Complete kit for quick degasification by suction of samples before analysis. It includes : a glass funnel filled with stainless steel wool, a 500 ml filter flask with stopper and a PP water-jet pump.

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PP water-jet pump, with sleeve nut R3/4", thread adapter R1/2", hose connection (olive 10/12 mm diam.) and vacuum connection (olive 6/9 mm diam.) removable with GL14 screwcap.

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Glass filtration unit. Complete kit with 300 ml funnel, sintered glass filter holder (porosity 10 µm), 1000 ml filter flask and aluminium clamp. Suction pump not included.

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3-position suction filtration unit. For membranes with diameter 47 or 50 mm. 300 ml stainless steel funnels, stainless steel filter holder, aluminium forceps and PTFE valves. Suction pump not included.

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Filter membranes, made of mixed cellulose esters. Suitable for microbiology, colony counting, sample prefiltration and clarification. Diameter 47 mm, box of 100 units.

	Porosité Capacity	Couleur / type Colour / type
650420	0,45 µm	Noire quadrillée / Black with grid
650410	0,45 µm	Blanche unie / Plain white
650411	0,65 µm	Blanche unie / Plain white
650412	0,8 µm	Blanche unie / Plain white
650413	1,2 µm	Blanche unie / Plain white

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Other sizes, porosities, colors : contact us

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Filter flask, with glass tube and olive.

	Capacité Capacity	Ø corps Body Ø	Ø col Neck Ø	Hauteur Height
617250	250 ml	85 mm	35 mm	145 mm
617500	500 ml	105 mm	35 mm	175 mm
617901	1000 ml	135 mm	45 mm	430 mm

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
LABORATORY GLASSWARE AND KITS

3.15. Pincettes et supports / Clamps and supports

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Conical rubber stoppers, with low taper ratio. They ensure a good airtightness and adhesion. Available without hole (P) or with hole (T).

	Øint	Hauteur
651109P/T	14/18 mm	20 mm
651110P/T	16/21 mm	25 mm
651111P/T	17/23 mm	30 mm
651112P/T	19/24 mm	25 mm
651113P/T	22/27 mm	25 mm
651114P/T	25/31 mm	30 mm
651115P/T	28/35 mm	35 mm
651116P/T	31/38 mm	35 mm
651117P/T	34/41 mm	35 mm
651118P/T	35/44 mm	40 mm
651120P/T	41/49 mm	35 mm
651121P/T	47/55 mm	40 mm
651122P/T	50/60 mm	45 mm
651124P	55/65 mm	50 mm
651126P	65/75 mm	55 mm
651127P	71/83 mm	60 mm
651128P	78/92 mm	50 mm
651133P	95/100 mm	55 mm
651137P	100/130 mm	85 mm

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Versilic hose, translucent, biologically inert. Complies with European Pharmacopeia chemical test requirements. Temperature of use : -50/+200°C. Possible peak at 220°C. Sold by the meter.

617600	Ø 1,5/3 mm	617602	Ø 5/7 mm
617605	Ø 3/7 mm	617604	Ø 6/9 mm
617601	Ø 4/6 mm	617603	Ø 8/12 mm

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Natural rubber hose. Possible sterilization in wet heat. Easy to fit on connectors thanks to its elasticity. Temperature of use : -35/+70°C. Sold by the meter.

617611	Ø 5/8 mm	617614	Ø 7/11 mm
617612	Ø 6/9 mm	617613	Ø 8/12 mm

Tuyau à vide. Gomme souple rouge, diamètre 8/18 mm. Vendu au mètre. Suction hose, made of red flexible rubber. Diameter 8/18 mm. Sold by the meter.

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Stainless steel rod with threaded base, with nut and washer.

600921	280 mm, Ø 8 mm
600922	380 mm, Ø 8 mm
600924	590 mm, Ø 10 mm

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Clamp with nut for tube and burette. Ni/Cr. Opening diameter 5 to 50 mm. Diameter of the rod < 18 mm.

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PP clamp with nut for 2 burettes. Diam. of the rod 10 to 14 mm.

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Complete kit for burette with one 170 x 110 mm base, a stainless steel rod 380 x Ø 8 mm and a clamp with nut.

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Clamp for stand without nut. Cork coated. Opening ≈ 40mm.

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Double nut for rods with maximum diameter 14 mm. Nickel.

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3-finger clamp with 2 tightening wheels. Chromium and zinc. Rubber covered jaws for a better grip.

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Open ring for stand. Ni/Cr. For funnels, flasks and iron gauze with ceramic gauze.

	Ø extérieur External Ø	Ø intérieur Internal Ø
600935	60 mm	46 mm
600930	100 mm	82 mm
600931	140 mm	125 mm

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- Accessoires de chauffage / Heating fittings

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
Silicon PTFE sealing ring. For SVL sliding joint.

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Silicon sealing ring. For SVL sliding joint.

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Rodavis screw assembly. Complete kit.

650162	19/26 mm
650161	29/32 mm
650160	24/29 mm

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PBT screwcap. GL range. Supplied with gasket.

650175	GL14	Diam.14
650176	GL25	Diam.25
650177	GL32	Diam.32

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PP connector for plastic hose. To assemble with GL14 screwcap

650180	Droit / Right
650181	Courbé / Curved

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Laboratory spatula. Robust and economical spatula, made of glass fiber filled nylon. Autoclavable.

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Plastic forceps, autoclavable. Length 118 mm.

3 594004 Pince inox pointe courbe 145 mm.

Stainless steel forceps, curved end. Length 145 mm.

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Glass weighing scoop, for standard use. With hollow stem to aid introduction of the powder into narrow neck vessels. The flat base ensures stability.

	Capacité Volume	Longueur Length	Matériau Material
594010	3 ml	70 mm	Verre / Glass
594011	6 ml	85 mm	Verre / Glass
594012	10 ml	100 mm	Verre / Glass

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7 427002 Goupillon standard / Standard brush

Standard drying rack, 55 positions. Made of PVC coated steel wire, with drip tray. Height 640 mm, width 360 mm.

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Hanging drying rack, 48 positions. Made of PVC coated steel wire. Height 500 mm, width 400 mm.

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Wall drying rack, 72 positions. Polystyrene. Equipped with 72 removable 100 mm pegs and a draining groove. Delivered with hanging hooks.

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Wall or hanging glass rack, chrome-plated. 5 rows for around 20 stemmed glasses in total.

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- Solution Glass Cleaner / Glass cleaner solution p.123

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From grape harvest, to vinification to ageing and bottling, the LQA can be used to analyze must, must under fermentation, and finished wine in less than 45 seconds.

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Built-in computer with a large 12-inch touch screen

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Tanks for automatic zeroing and cleaning cycles

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Analytical capacity of 60 samples per hour

Multiparamètres Multiparameters

Alcool, sucres, acidité totale, acide malique, pH, densité...

Alcohol, sugar, total acidity, malic acid, pH, density...

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Must, must under fermentation, or finished wine

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Technical sheet and demonstration video on our website

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Magnetic mini-stirrer. Variable speed, from 100 to 1000 rpm. For small volumes (maximum 1 litre). ABS tray, 120x120 mm.

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Magnetic stirrer. Analogic adjustment of speed, from 100 to 2000 rpm. For stirring up to 15 litres. Stainless steel tray, 150x150 mm.

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Heating magnetic stirrer. White glass ceramic top with excellent chemical resistance, easy to clean. Cast aluminium body shaped for stability.

- Analogic adjustment of speed (maximum 2000 rpm) and temperature (maximum 450°C). Powerful stirring up to 15 litres.
- Dimensions of tray / heated area : 150 x 150 / 120 x 120 mm
- Fast heating and accurate temperature stabilization
- «HOT» warning light when the plate temperature is above 50°C, even if the device is off.

Fitting for retort rod (not included).

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Mini vortex mixer. Variable speed 2500-4500 rpm with simple push to operate control. Perfect for mixing in cuvettes for spectrophotometry and tubes of maximum diameter 30 mm. Weight 550 g.

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Stirring bar retriever. PTFE-coated magnetic stem length 330 mm.

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Bench autoclave. Low footspace.

Permanent indication of the cycle steps.

Water level monitoring of the reservoir and sterilization chamber. Water tank 4,2 litres.

Safety : automatic control of the water level in the tank before the cycle starts ; over-heating, safety valve and overpressure switch.

	Volume de la chambre Chamber volume	Puissance Power	Temp. Temp.	Operation Modes
606501	16 litres (Ø230 x 410 mm)	1,9 kW	121 ou/ou 135°C	8 programmes pré-installés 8 pre- registered programmes
606502	24 litres (Ø260 x 450 mm)	2,4 kW		

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BagMixer 400CC. Mixer for samples of 400 ml. Particularly adapted to grape blending, for measurements on must before fermentation, thanks to an efficient work. Glass door, click-on paddles for an easier cleaning.

Volume 50 to 400 ml. Adjustable speed and time.

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Compact and economical laboratory balance. Excellent value for money. Large LCD. «TARE» function. PVC weighing plate. 9V battery.

	Portée Weighing range	Lecture Readout	Plateau Weighing plate
601707	500 g	0,1 g	Ø 150 mm
601704	2200 g	1 g	Ø 150 mm
601708	5200 g	1 g	Ø 150 mm

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Precision balance. Stainless steel weighing plate. Delivered with power supply adapter, control weight and protective cover. RS232 output for printer or PC. Integrated programmes: parts counting, percentage determination.

	Portée Weighing range	Lecture Readout	Plateau Weighing plate
601700	200 g	0,01 g	Ø 105 mm
601702	400 g	0,1 g	Ø 105 mm
601705	2000 g	0,1 g	130 x 130 mm
601706	6000 g	1 g	150 x 170 mm

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Compact analytical balance. Modern all glass draft wind shield for optimal visibility of the weighing object. GLP/ISO record keeping of weight values with date, time and serial number. Higher weighing speed by switching off the last digit. Totalising of weights (manual/automatic). Possible homologation in legal metrology.

Alibi memory: paperless archiving of up to 100,000 weighing results. Internal memory for up to 999 weighing results, 1000 items or recipe ingredients, 100 container weights, 100 users.

Dispensing mode: by pressing one key you can set all relevant parameters for the dispensing process.

Weighing range: 160 g; readout 0,1 mg
Adapter 230V. RS232 and USB output.

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
Other balances : contact us.

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LABORATORY EQUIPMENT

4.4. Balances (suite) / Balances (continued)

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Scale with large platform. Simple and convenient 4-key operation. Wall mount for display device. Stainless steel platform, 55x55x6 cm. Power supply 230 V or batteries.
Measurement range: 0-300 kg; readout 0,1 kg.

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Floor scale. Weighing bridge made of anti-slip corrugated steel, 4 silicone-coated steel load cells. Easy levelling. Large display (25 mm high figures). Universal power adapter. Optional ascending ramp or stable pit frame. Possible homologation in legal metrology.

	Portée Range	Lecture Readout	Charge min. Min. load	Plateau Weighing plate
601750	600 kg	0,2 kg	4 kg	120x150x11 cm
601751	1500 kg	0,5 kg	10 kg	120x100x11 cm

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Hanging scale. Perfect to weigh bags and buckets. Hold feature. Large display (12 mm high figures). 3 AA 1,5 V included batteries.

	Portée Weighing range	Lecture Readout
601760	50 kg	20 g
601762	100 kg	200 g

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DKD calibration (German calibration service). Valid internationally, it holds traceability to the national standard and determines the deviation from true, actual measure value.

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Compact ultrasonic bath. Capacity 0,75l. Digital control of timer : 90-180-280-380-480 sec. Bath opening (L x l x P) : 15 x 13 x 5 cm. Supplied with basket and lid.

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Ultrasonic bath. Indispensable at the laboratory to degas samples, to clean light pieces, etc. Rugged industrial transducers 40 kHz, a frequency that eliminates standing waves and creates an homogeneous cavitation in the tank.

Accessible control tower located behind and above the bath for safer operation.

Permanent operation/supply control, to maintain ultrasonic power, depending on the liquid level and temperature changes.

Drain valve with hose nipple and 1 m hose.

Plastic body, resistant to chemical products. Plastic tank cover included.

Mechanical timer, which can be set for up to 60 minutes and has a setting for continuous cleaning.

	Capacité Volume	Dimensions cuve (LxPxH) Tank dimensions (WxDxH)
609000	5,7 litres	29 x 15 x 15 cm
609001	9,5 litres	29 x 24 x 15 cm

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Additionally to the characteristics detailed above, this model is equipped with a 99-minute digital timer and continuous ultrasonic operation. High/low power control for more or less delicate applications.

	Capacité Volume	Dimensions cuve (LxPxH) Tank dimensions (WxDxH)
609005	5,7 litres	29 x 15 x 15 cm
609006	9,5 litres	29 x 24 x 15 cm

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PURA10 water bath. Equipped with a microprocessor maintaining a temperature stability of $\pm 0,15^{\circ}\text{C}$, this equipment offers many application possibilities. LED display and keyboard with splash-water protection. Easy-to-use. Lid not included.

Working temperature range : Ambient to $+99,9^{\circ}\text{C}$.

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Transparent lift-up lid for PURA10 water bath, out of plexiglas and polycarbonate.

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SW22 shaking water bath. Equipped with a microprocessor for a $\pm 0,2^{\circ}\text{C}$ temperature control, this equipment is conceived for unpervised operation : electronic timer, high/low temperature limits, automatic shut-down in case of low level. Splash-proof LED display and keyboard. Cover not included.

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Transparent lift-up Makrolon cover, for SW22 shaking water bath.

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Pompe refoulante : 6 l/min, 0,1 bar.

Open Heating Bath Circulator C BT9. The circulator is mounted onto a bridge and can be easily removed. Transparent bath tank. Bright white display, easy to read. Precise temperature regulation, with stability $\pm 0,03^{\circ}\text{C}$. Very quiet and easy to use.

Working temperature : ambient to $+100^{\circ}\text{C}$.

Bath opening (L x W x D) : 23 x 15 x 15 cm (Volume 6-9 litres)

Waste pump : 6 l/min, 0,1 bar.

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LABORATORY EQUIPMENT

4.6. Bains thermostatés et cryostats (suite) / Refrigerated - heating baths (continued)

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Refrigerated heating circulators to heat or cool external systems (refractometers, photometers, viscosimeters) and an internal bath tank made of stainless steel for small samples. High heating/cooling capacity to quickly reach the setpoint. Convenient cleaning of dust by removable venting grid. All models provide a drain port for conveniently draining the bath fluid.

- Bright white display temperature display, resolution 0.1 °C
- Lid included
- Temperature stability $\pm 0,03^\circ\text{C}$
- Waste pump : 15 l/min, 0,35 bar.

	Temp. de travail Working temp.	Capacité bain (LxPxH - volume) Bath capacity (WxDxH - volume)
608040	-20...+100°C	13 x 15 x 15 cm / 3...4 litres
608050	-38...+100°C	35 x 41 x 30 cm / 48...56 litres

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For all water baths, refrigerated and heating circulators with baths, many accessories are available (contact us) :

- Bath tanks to connect to circulators
- Tubes holder (different sizes)
- Adjustable platforms
- Bath fluids

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Dry heating block for 10 cuvettes. For the incubation of mixtures in cuvettes with optical path 10 mm before measuring the optical density. Temperature range : +25°C...+42°C ($\pm 0,1^\circ\text{C}$)

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Weight ring, for the stabilization of Erlenmeyer flask in an ultrasonic or thermostatic bath.

	Poids Weight	Ø intérieur Internal Ø	Pour fioles For flasks
606015	220 g	48 mm	125-500 ml
606016	450 g	56 mm	250-1000 ml
606017	680 g	61 mm	500-2000 ml
606018	910 g	74 mm	1 - 4 litres

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Mini-centrifuge Frontier 5306. Small and powerful device for simple and quick liquid separation, this centrifuge provides safe and easy operation for everyday applications. Brushless quiet motor, auto shut-off safety switches. Supplied with an 8 place 1.5/2 ml tube rotor as well as a 4 place PCR 8-strip tube rotor. Maximum speed 6000 rpm.

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2 ml Clearlock microtubes (x1000), with stopper.

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Bench centrifuge Frontier 5706. High capacity in a compact structure to optimize the space on laboratory benchtop. Quality raw materials. Digital bright display and intuitive interface.

Programmable centrifugation time from 10 sec. to 99:59 minutes or continuous. Automatic locking system of the lid (motorized). Brushless quiet motor. Imbalance detection system with automatic stop.

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Fixed angle rotor 6 x 50 ml. For tube with diameter 31 mm.

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Glass alcohol lamp, 100 ml. With safe system, for wick Ø 4 mm.

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Flask holder, tripod, made of galvanized steel. Height approx. 21 cm.

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Iron wire gauze, made of galvanized iron wire mesh, with ceramic center.

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Gas burner with stopcock and pilot flame. Tube diameter 11 mm. For propane/butane gas.

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High flow burner Ø 40 mm. Ideal to heat big capacity water tanks (1000 ml).

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Gas burner holder with nut. It allows to adjust the height of the burner ref. 640123, fixing it on a rod with maximum diameter 12 mm.

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Metal-cased heating mantle. Plastic coated, chemical-resistant housing with built-in power-on and heating zone switch. Flexible glass yarn heating element. 2 heating zones (min/max). Built-in power regulator. Max. heating temperature 450°C. 1.5 m power supply (earthed) cable. Nominal voltage 230 V.

	Capacité Capacity	Puissance Heating power
606020	250 ml	150 W
606021	500 ml	200 W
606022	1000 ml	300 W
606023	2000 ml	500 W

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6-place heating mantle. Same structure as the above model, only one heating zone per position. Built-in power regulator with analogic adjustment at each position. Dimensions : 1150 x 150 x 290 mm. Power supply : 230 V.

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Clamp for rod with diameter up to 12 mm. Appropriate for heating mantles ref. 60602X and 60605X.

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Stainless steel elevator. Steady and stable structure. 250 x 250 mm tray, adjustable height from 65 to 270 mm. Maximum weight 3 kg.

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3-place heating mantle, or heating-mantles with other capacities, with magnetic stirrer: contact us.

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Dimensions (H x L x P), poids : 50 x 15 x 45 cm, 10 kg

The BASIC PH4 water distillation unit is a continuous distillator. Easy to run, it is equipped with a safety thermostat. Production : 4 litres per hour thanks to a 3 kW chromium plated heater. The horizontal boiler is made of high quality borosilicate glass. CE certification and marking.

*Conductivity (raw water) : 3.0 - 4.0 µs/cm
Conductivity (treated water) : 1.5 - 2.0 µs/cm
PH : 5.5 - 6.0*

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Incubator CLN. Natural air convection. Stainless steel interior, plain door. Controlled by microprocessor. Digital display. Safety thermostat. Alarm in case of faulty temperature probe. Auto-diagnosis function. Control certificate at 37°C. RS232 output. Temperature range: 5°C above ambient temperature...+100°C Uniformity : ± 0,5°C (except for 56-liter volume : ±0,2°C)

	Volume int. Int. volume	Finition Finish	Etagères (standard / maxi) Shelves (standard / maxi)
606101	15 litres	inox / stainless steel	1 / 2
606102	15 litres	epoxy blanc white epoxy coated	1 / 2
606111	32 litres	inox / stainless steel	1 / 3
606112	32 litres	epoxy blanc white epoxy coated	1 / 3
606121	56 litres	inox / stainless steel	2 / 5
606122	56 litres	epoxy blanc white epoxy coated	2 / 5

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Incubators with higher capacity, with refrigeration : contact us.

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Pocket magnifier with folding handle. Glass lens Ø40 mm. Chromium plated brass frame and handle. Magnifying factor: 4x

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Pocket lens made of organic glass. Black plastic case. Dimensions: 53 x 53 mm. Magnifying factor: 3,5x.

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Mini-microscope. This pocket instrument allows to visualize with precision the vine parasites thanks to an integrated lighting and to an appropriate magnifying factor. Magnifying factor : 30x.

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Stereomicroscope 20x/40x/80x, reflection/transmission type. Equipped with a 360° rotating head with 2x and 4x objectives. WF eyepieces 10x and 20x. LED illumination. Includes a dustcover, a pair of rubber eyecups, power supply adapter, glass plate and plastic plate (white/black).

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Researcher binocular microscope 40x-1000x. Transmission microscope adapted to school, study and laboratories. Nosepiece with DIN objectives 4x, 10x, 40x and 100x oil (achromatic). Eyepieces pair : WF 10x (Ø23mm). 20W halogen illumination with dimmer. Abbe condenser (N.A. 1.25). Includes a dustcover, immersion oil and power cable. Trinocular model also available.

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Bioscience trinocular microscope 40x-1000x achromatic. Transmission microscope adapted to study and science. The ergonomic design facilitates its use. The trinocular head makes possible the simultaneous connection with the optional Mikrocam camera, for a permanent and simultaneous observation on screen. Nosepiece with DIN lenses 4x, 10x, 40x and 100x oil (achromatic). Eyepieces pair : WF 10x (Ø23mm). Coaxial mechanical desk with nonius scale. Coarse and precise focusing. 30W halogen illumination with dimmer, Köhler type. Abbe condenser (N.A. 1.25). Includes a dustcover, immersion oil and power cable.

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WF 20x eyepiece (Ø23mm), for Bioscience microscope.

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Mikrocam 1,3 MP 1280 x 1024. High quality of images. Possibility to make video sequences and single images. PC connection through USB 2.0 port. Professional software with functions of image capture and measurement. Adjustable white balance and exposure control, manually or automatically. Supplied with 0.5x reducing lens and adapters.

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Wide range of available accessories for microscopes : contact us.

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
This furniture range offers a wide choice of dimensions, worktop materials, end pieces, fittings, etc. It has been conceived to be set up easily. The parts are delivered in spare parts but can be mounted upon request.

The elements hereafter are only a selection.

For an on-site or distance study, or to obtain information about the whole range (other coating/materials, other kind of furniture): contact us.

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16 m² laboratory approximately


Suggested furniture (blue parts) :

- Dry benches with stratified work surface and top shelf on support, on metal leg
- Wet bench with stratified work surface and top shelf on support, on metal leg, equipped with a laboratory PP sink and epoxy-coated fittings
- Shelf under bench
- Stoneware sink unit with epoxy-coated fittings on furniture with top shelf on support
- Cabled 4 and 3-plug blocks
- Wall-mounted cupboard with solid doors
- Mobile cabinet with 4 drawers
- Mobile cabinet with 2 doors and a shelf
- Mobile cabinet with 2 doors and 1 drawer
- Laboratory chair on wheels
- Support for computer

MOBILIER DE LABORATOIRE

LABORATORY FURNITURE

5.1. Suggestions d'aménagement (suite) / Layout suggestions (continued)

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40 m² laboratory approximately


Suggested furniture (blue parts) :

- Dry benches with stratified work surface and top shelf on support, on metal leg
- Wet benches with stratified work surface and top shelf on support, on metal leg, equipped with a laboratory PP sink and epoxy-coated fittings
- Stoneware sink unit with epoxy-coated fittings on furniture with top shelf on support
- Central benches with technical box, stratified work surface
- Cabled 3-plug blocks
- Wall-mounted cupboard with solid doors
- Mobile cabinet with 4 drawers
- Mobile cabinet with 2 doors and 1 drawer

MOBILIER DE LABORATOIRE

LABORATORY FURNITURE


5.1. Suggestions d'aménagement (suite) / Layout suggestions (continued)

MOBILIER DE LABORATOIRE

LABORATORY FURNITURE

5.1. Suggestions d'aménagement (suite) / Layout suggestions (continued)

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45 m² laboratory approximately including :

- An area for chemical analysis with office section, 40 m²
- An area for microbiology, 6 m²

Suggested furniture (blue parts) :

- Dry benches with stratified work surface and top shelf on support, on metal leg
- Wet benches with stratified work surface and top shelf on support, on metal leg, equipped with a laboratory PP sink and epoxy-coated fittings
- Stoneware sink unit with epoxy-coated fittings on furniture with top shelf on support
- Central benches with technical box, stratified work surface
- Safety microbiological unit with laminar flow equipped with a UV lamp
- Cabled 3-plug blocks
- Wall-mounted cupboard with glass doors
- Mobile cabinet with 4 drawers
- Mobile cabinet with 2 doors and 1 drawer
- Built-in cabinet under bench with sliding doors
- 2-body cabinet (bottom with solid doors, top with glass doors)

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100 m² laboratory approximately including :

- An area for chemical analysis with office section, 50 m²
- An area for microbiology, 15 m²
- An area for quality control, 15 m²
- An area for washing, 6 m²
- An area for sample storage, 5 m²

Suggested furniture (blue parts) :

- Dry benches with stratified work surface and top shelf on support, on metal leg
- Wet benches with stratified work surface and top shelf on support, on metal leg, equipped with a laboratory PP sink and epoxy-coated fittings
- Stoneware sink units with epoxy-coated fittings on furniture with top shelf on support
- Central benches with technical box, stratified work surface
- Laboratory fume hood equipped with a fan system
- Safety microbiological unit with laminar flow equipped with a UV lamp
- Anti-vibration table for analytical balance with top shelf on support
- Cabled 3-plug blocks
- Wall-mounted cupboard with glass doors
- Mobile cabinet with 4 drawers
- Mobile cabinet with 2 doors and 1 drawer
- Mobile cabinet with 1 drawer and 1 file drawer
- Cupboard with glass doors on both faces on technical box
- Built-in cabinet under bench with sliding doors
- 2-body cabinet (bottom with solid doors, top with glass doors)
- High cabinet with 2 vertical drawers
- Standard cabinet with 2 doors 4 shelves

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Multirisk safety cabinet, for the storage of non-flammable chemical products (acids, bases and solvents with low inflammability degree). Completely manufactured from 1,5 mm thickness steel sheet, cold pressed, acid resistant epoxy coated. One-piece structure with screw-jacks for level adjustment.

- Separation for "product category" in 3 compartments.
- Locking doors.
- Doors opening angle 110° to remove shelves without leaning them.
- Ring Ø125 mm for external connection or extraction system
- PP external hinges, corrosion resistant.
- Compartments equipped with basin form shelf, height adjustable, with welded corners.
- Safety signs complying with ISO 3864 norm.

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Accessories for multirisk safety cabinet.

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COMBISTORAGE acids/bases safety cabinet. Underbench version of multirisk safety cabinet. 2 compartments, 2 doors, 990 x 520 x 620 mm.

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Safety cabinet for flammable products. Firewall structure. Doors built with sleeves of various materials, to ensure a good resistance to fire. Resistant to high temperatures (800 °C).

- Automatic door closure in case the ambient temperature exceeds 50 °C ; Lock with key and door blockage ; Anti-sparking hinges.
- 30 mm thermo-expanding gasket (DIN 1402).
- Certified valves close the ventilation holes when the temperature exceeds 70 °C.
- Air extraction duct collar 100 mm diam. for a connection to outside or to a filtering aeration system.
- Grounding connection.
- Feet adjustment level.
- 3 epoxy coated shelves on stainless steel supports, height adjustable ; bottom spill pallet with draining grid.

- 853101** Armoire 1 porte, 3 étagères, bac de rétention et caillebotis. Dimensions : 595 x 600 x 1950 mm
Cabinet with 1 door, 3 shelves, retention container and spill pallet. Dimensions : 595 x 600 x 1950 mm
- 853102** Armoire 2 portes, 3 étagères, bac de rétention et caillebotis. Dimensions : 1200 x 600 x 1950 mm
Cabinet with 2 doors, 3 shelves, retention container and spill pallet. Dimensions : 1200 x 600 x 1950 mm

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Accessories for safety cabinet for flammable products.

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COMBISTORAGE safety cabinet for flammable products. Underbench version of safety cabinet for flammable products, 1 drawer. Dimensions : 990 x 520 x 620 mm.

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Vertical laminar flow cabinet. Conceived for microbiological and cell culture applications, it protects products or growth media (non-pathogen) in sterile environment. Epoxy-coated steel frame. Perforated stainless steel (AISI 304) worktop, 2 panels, removable for cleaning. Temperate glass front window, possibility to lock in the fully open position. Thermally protected motor fan, with manual regulation on control panel with "soft touch" membrane keyboard. Noise level: <58 dB(A). Absolute HEPA filter H14 class and pre-filters. Warning visual on display of exceeding the maximum period for the replacement of filters. Digital counter of working time. Illumination work area > 700 lux (13W fluorescent white lamp). External dimensions (L x W x H) : 800 x 580 x 940 mm. Power supply : 230V - 50 Hz.

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Fume hood, EN"1" fire resistance class. Designed for chemistry laboratories, to protect operators from breathing harmful vapors of chemical products. Wood water/fireproof high-pressure laminate panels and working area. Except for the metal lateral frame, acid-resistant epoxy coated. Plain sides. Back panel made of PP, removable for a total cleaning of the fume hood. Low footprint on the bench. The side airfoil service panels are used for electric utilities on the left, and fluids on the right (water, gas, etc.). Mobile front sash. Made of Securit® glass 6 mm. Stop at 400 mm of the worktop. Break system on the sides. Air flow optimization towards the centre of the fume hood to prevent vapors from going out.

Extraction system to be ordered separately. Extraction outlet diam. 250 mm. Total depth : 870 mm ; height 1600 mm. Minimum height for installation : 2700 mm. Fluorescent lamps 400 Lux - IP65. Complies with normes EN14175 and XPX15206. Classification : NFX15210 class B.

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Fume hood with 1 front sash to be set up on existing bench.

- 851012** Largeur 1200 mm / *Width 1200 mm*
- 851015** Largeur 1500 mm / *Width 1500 mm*
- 851018** Largeur 1800 mm / *Width 1800 mm*

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Floorstanding fume hood with 1 front sash. Metal support frame «C» shape, height 950 mm. Modular legs with epoxy-coated steel tubes. Adjustable screw-jack (minimum height 87 cm). Glass worktop. Optional back and side panels.

- 851112** Largeur 1200 mm / *Width 1200 mm*
- 851115** Largeur 1500 mm / *Width 1500 mm*
- 851118** Largeur 1800 mm / *Width 1800 mm*

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Test strips (x100). Colorimetric reagents. Results obtained by comparison with a colorimetric scale.

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Free and total chlorine portable photometer, wide range. LED light source, filter at 525 nm. Calibration control function. Intuitive use and measurements very easy to run. Delivered with 2 cuvettes and 3 1,5V AA batteries. Reagents not included.
Measurement range : 0,06 to 10,00 mg/l. Accuracy : ± 0.03 mg/l $\pm 3\%$ of reading ; excluding dilution error.

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Hardness (total) portable photometer. LED light source, filter at 466 nm. Calibration control function. Intuitive use and measurements very easy to run. Delivered with 2 cuvettes and 3 1,5V AA batteries. Reagents not included.

Measurement range :
Narrow range : 0 to 250 mg/L ; ± 5 mg/L ± 4 % of reading
Medium range : 200 to 500 mg/L ; ± 7 mg/L ± 3 % of reading
Wide range : 400 to 750 mg/L ; ± 10 mg/L ± 2 % of reading

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The implementation of a hygiene plan, incorporating ATPmetry A3 as a management tool, is now recommended for all wineries, not only to ensure better wine quality but also to streamline cleaning processes, which can lead to reduce water and chemical consumption.

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Lumitester Smart ATP-meter. This portable and light instrument gives practical information about the degree of hygiene of surfaces. The sample is taken by wiping a pad on a bounded surface, then this pad is inserted into a tube where it reacts with enzymes. This tube has to be inserted into the instrument to measure bioluminescence produced by this reaction.

Hygiene monitoring in three easy steps : Swab the sample, measure and analyse. Results within 10 seconds. Centralized multi-site data with cloud data management. Visualization and permanent control of data thanks to the application. Measurement in RLU – Relative Light Units. 2000 storage spaces. Delivered with 2 AA batteries, cleaning brush and USB cable.

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Swabs for Lumitester Smart ATP-meter with A3 reagents. Larger detection of residues and micro-organisms thanks to the conversion of AMP and ADP into ATP. The use of these tests including a triple detection (A3) is an asset to take into account all kinds of living cells and avoid false negatives. Supplied by 100 pieces, in 5 packs of 20 tubes containing enzymes and a swab. To be kept cold (2-8°C) for one year.

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Protection against pollution FFABE1P3 :

A1: protection against gas and organic vapors (boiling point > 65°C). Ex.: painting solvents, adhesives, pesticides.

B1: protection against gas and inorganic vapors (except carbon monoxide). Ex.: acetone, butane.

E1: protection against acids, sulphur dioxide (SO₂)

P3: protection against toxic dust and particulates, smoke, fogs (solid particles or aqueous aerosol). Ex.: copper, welding fumes, sodium hydroxide, calcium oxide, fibers, micro organisms, biochemical substances.

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Temporary half-mask (1-month use). FFABE1P3 protection, high performance. Integrated filters with flat shape for wider field of vision. Front piece made of thermoplastic material, hypoallergenic. Low resistance to respiration.

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Preassembled filters : combination of a P3 dust filter and a A1B1E1K1 gas filter. For Easylock filter-holder. Box of 6 pieces.

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Latex, ambidextrous single-use gloves. Non sterile, powder-free, textured fingertips. Rolled cuffs minimise the risk of tearing. Suitable for food handling. Color : natural. Box of 100.

- 199700** Taille 6/7 (S) / Small size (6/7)
- 199701** Taille 7/8 (M) / Medium size (7/8)
- 199702** Taille 8/9 (L) / Large size (8/9)

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Durable apron made from strong PVC. With nylon neck band. High resistance to oils, greases and certain chemicals. Dimensions : 84 × 112 cm.

6.3.3. Premiers secours / First aid equipment

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Eye wash station, containing 2 solutions in flask to provide a hygienic and quick eye wash in case of contamination or irritation: «pH Neutral», a phosphate buffer solution for to neutralize acids and alkalis ; «Eye Wash», a sodium chloride solution to rinse dust and dirt. Ergonomic eye cup on the flask, for a hygienic and regular application.

	Forme Type	pH Neutral pH Neutral	Lave-œil Eye wash
199750	Station murale Wall station	200 ml	500 ml
199751	Coffret mural étanche Waterproof wall case	200 ml	500 ml
199752	Station murale Wall station	500 ml	1000 ml
199753	Coffret mural étanche Waterproof wall case	500 ml	1000 ml

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Closed glass wine thief with rubber handle. Round end with a hole to be stoppered with the thumb. Equipped with a rubber handle, food grade.

Droit / Straight	Coudé / Curved	Longeur / Length
501025D	-	25 cm
501035D	501035C	35 cm
501045D	501045C	45 cm
501055D	501055C	55 cm
501065D	501065C	65 cm
501075D	501075C	75 cm
501095D	501095C	95 cm

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Personalized marking (one or more colors) : contact us.

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Open wine thief with rubber handle. Open end with a food grade removable rubber stopper with a hole to be stoppered with the thumb. Possibility to brush the inner surface with a bottle brush. Equipped with rubber handle, food grade.

Droit / Straight	Coudé / Curved	Longeur / Length
501125D	-	25 cm
501135D	501135C	35 cm
501145D	501145C	45 cm
501155D	501155C	55 cm
501165D	501165C	65 cm
501175D	501175C	75 cm
501195D	501195C	95 cm

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SAMPLING TOOLS AND CONTAINERS

7.1. Tête-vin (suite) / Wine thieves (continued)

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Straight wine thief with 4 holes, equipped with a rubber handle, food grade.

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Closed wine thief with glass handle(s). Round end with a hole to be stoppered with the thumb.

Droit / Straight	Coudé / Curved	Longeur / Length
504045	503045	45 cm
504055	503055	55 cm

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Glass tube for spirit sampling (topette).

507000	50 ml environ, diamètre 30 mm 50 ml approximately, diameter 30 mm
507002	100 ml environ, diamètre 40 mm 100 ml approximately, diameter 40 mm
507001	170 ml environ, diamètre 40 mm 170 ml approximately, diameter 40 mm

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Sample press. Adapted to fleshy fruit without stone. Equipped with a handle, aluminium structure, stainless steel cage and a drainpipe.

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Stainless steel ruler for barrel. Double graduation. With hanging ring.

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Polypropylene jug, conical shape. Autoclavable at 121°C. Ergonomically designed handle.

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Aspirator bottle with spigot, made of high density polyethylene. Heavy duty material (2mm wall), conceived to be used for a long time. It includes a robust spigot to ensure a perfect leakproof performance. A few drops of vaseline shall be put on the thread of the spigot before assembly. 1 or 2 handles according to capacity.

	Capacité Volume	Diam.col Neck diam.	Diamètre Diameter	Hauteur Height
803001	5 l	46 mm	163 mm	335 mm
803002	10 l	55 mm	210 mm	414 mm
803003	25 l	55 mm	279 mm	525 mm
803004	50 l	55 mm	366 mm	603 mm

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Tank with spigot. Made of high density polyethylene, food grade. Specifically designed to rationalize laboratory organization. 1 handle.

	Capacité Volume	Dimensions Dimensions
803011	10 l	380 x 300 x 160 mm
803012	20 l	470 x 370 x 180 mm

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Ribbed cap PP 28. For syrup bottle. Self seal.

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Glass flask, borosilicate glass. Ergonomic shape. GL45 opening.

	Capacité Volume	Diamètre Diameter	Hauteur Height
801020	125 ml	55 mm	124 mm
801021	250 ml	66 mm	158 mm
801022	500 ml	78 mm	193 mm
801023	1000 ml	93 mm	253 mm

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Invio 25 cap. Black, made of PE, with tamper-proof ring. For square bottle with handle.

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Invio 25 cap with dropper. White, made of PELD. For square bottle with handle.

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100ml dropping bottle, PE-LD, GL18 dropping insert and locking cap.

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PP cup. With blue screw cap. Capacity 150 ml.

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SAMPLING TOOLS AND CONTAINERS

7.6. Flaconnage (suite) / Bottles (continued)

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Wash bottle in polyethylene. Tube with dispensing tip which gives ultra-fine stream, removable to increase flow.

	Capacité Volume	Diamètre Diameter	Hauteur Height
601050	50 ml	35 mm	96 mm
601100	100 ml	43 mm	115 mm
601250	250 ml	60 mm	140 mm
601500	500 ml	75 mm	180 mm
601901	1000 ml	95 mm	220 mm

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Printed wash bottle, «Ethanol». Wide mouth with a valve that prevents liquid drips. Capacity 500 ml.

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Printed wash bottle, «Distilled water». Wide mouth. Capacity 500 ml.

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Our standard solutions and reagents are prepared and packed according to protocols described in the Quality Assurance Management System of the company.

The standard solutions are systematically controlled and adjusted before conditioning, using stabilized standard solutions.

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REAGENTS AND SOLUTIONS

8.1. Solutions titrées (suite) / Standard solutions (continued)

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Acidimetric solution for total acidity.

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Acidimetric solution for volatile acidity.

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Acetimetric solution for acetic acid.

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Blue sodium hydroxide. With bromothymol blue. Concentration N/10. Short validity (1 month, 15 days after opening).

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916000B 500 ml Fehling B

Iode / Iodine.

Concentration N Concentration N/10 Concentration N/50

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Validity short (1 month, 15 days after opening).

Mixed solution for chromatography 250 ml. Short validity (1 month, 15 days after opening).

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931250 250 ml 931500 500 ml 931901 1000 ml

Soude (hydroxyde de sodium) / Sodium hydroxide.

Concentration 1,125N Concentration 2N Concentration 4N

934250 250 ml 908250 250 ml 909500 500 ml
934500 500 ml 908500 500 ml 909901 1000 ml
934901 1000 ml 908901 1000 ml

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Colored reagent for Frantz Paul method, 500 ml.

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Milk of lime. Validity 2 years, 1 year after opening.

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Antifoam. Validity 1 year, 3-4 months after opening.

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Glass cleaner, ready-mixed alkaline cleaner for laboratory glassware.

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Titrvin offers a complete range of reference wines to calibrate apparatuses; and to assess and control methods of measurement or analysis. Complying with the ISO Guide 34, the batches are prepared by a COFRAC accredited laboratory for the activities of Reference Material Producer.

In order to guarantee the «wine matrix» effect, each reference wine is exclusively prepared by mixing wines and musts from various origins.

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The Titrvin AA series offers a range of ordered values of the 9 most common parameters to calibrate the response patterns of automatic equipments. Packaged in box of 20 glass vials of 10 ml.

100511	AA1	100514	AA4
100512	AA2	100510	AA5
100513	AA3		

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The Titrvin SW series offers a range of ordered values of the 9 most common parameters to calibrate the response patterns of automatic equipments. Increasing alcohol content and decreasing sugars content, designing a matrix of fermenting musts or wines with residual sugars. Packaged in box of 20 glass vials of 10 ml.

100521	SW10
100522	SW20
100523	SW30

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BTA Titrvin. Control standard for around 15 parameters. Used every 10 samples of wine analysed, it helps to control stability of measures on automatic analyzers. Box of 20 vials.

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BTB Titrvin. Control standard for manual analysis, with around 15 determined parameters. To be used periodically to validate the measurement and analysis methods run at the laboratory. Box of 6 bottles 250 ml.

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Ethycal, hydro-alcoholic solution with certified value of density. 250 ml flask, delivered with calibration certificate linked to a COFRAC certificate, giving the density and the calculated AbV. Approx. 12% AbV.

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Standard wine. Recommended for the calibration of electric ebulliometers. AbV determined by OIV official method. Certificate of analysis provided upon request.

100552	250 ml
100550	500 ml

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Malic acid, 3 standard samples for the determination of malic acid by chromatography on Whatman paper. 10 ml flask.

943601	1 g/l
943602	2 g/l
943603	3 g/l

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Calibration solution for refractometer. 60 ml flask.

196750	10% Brix
196751	15% Brix
196752	20% Brix

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Other concentrations up to 60% Brix : contact us.

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Practical alcohometric guide, containing tables with :

- Corrected alcohol content in % AbV according to the AbV reading and temperature. Range : -10 to +40°C (1/2°C) and 0 to 100% AbV (1/10 % AbV)
- Pure alcohol volume at 20°C in 100 dm³ of a hydro-alcoholic solution
- Content of methyl alcohol (% volume)
- Content of propyl alcohol
- Densities
- Calculations to increase/reduce the alcohol content of spirits.

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Alcohometric tables expressed in degrees. Alcohol content in % AbV according to the AbV reading and temperature. Range : 0 to +40°C (1°C) and 0 to 100% AbV (1% AbV).

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Gauging scale. Capacity of barrels, any shape.

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Inventory triple booklet. Capacity and volumes contained in containers. Extended range.

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Correspondence table Baumé / Specific gravity.

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
p.131

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BOOKS AND SOFTWARES

9.2. Livres techniques / Technical books

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BOOKS AND SOFTWARES

9.2. Livres techniques (suite) / Technical books (continued)

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LABOX, digital cellar assistant. Applications for daily process support and calculations: volume and AbV conversions, barrel or tank volume calculations, reductions, blends...

- Alcoholometry: AbV and volume conversions, barrel volume calculations, pump-overs and reductions, fortification, SO₂, energy value.
- Inventory: stock file with or without wet and AbV measurements; Pure alcohol and losses calculated with a single click.
- Digitank: gauging tables for tanker trucks, barrels, or tanks. Pure alcohol to remove, gauge height to reach, regular volume readings.
- Optimix: proportions of different products to blend to achieve a solution that meets cost, analysis, or volume objectives.

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Metrology is an essential part in product quality. It is a pillar in the same way as normalization and certification. As a matter of fact, accuracy of measuring instruments is fundamental to control manufacturing processes and assure the conformity of analyses.

Any measuring instrument is subjected to modifications due to wear and tear or clogging. Its metrological performances evolve, in a non-predictive way. Maintenance and regular control of measuring instruments allow to assure good results, to respect recommendation of ISO 9001:V2015 norm and to control procedures.

Laboratoires Dujardin-Salleron have set up the metrological department with the company E.R.IN (Study and Research in Industry), with the purpose of calibration, control, verification and maintenance of instruments. The provided services comply with norms in force and lead to issue a certificate that guarantees traceability to international system of units.

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REFRACTOMETRY : FROM 0 TO 72,8% BRIX

Type : handheld refractometer, optical or digital
Recommended frequency : on an annual basis, or after a long period without use
Self-control : distilled water or standard solutions (cf.p.125)
Metrological services : calibration and adjustment in 2 points
Method : use of a standard block or calibration solutions

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THERMOMETRY : -20 to +150°C

Type : expansion (stem or enclosed type) or digital
Recommended frequency : on an annual basis for electronic thermometers, and according to internal norms of the user institution for expansion thermometers
Other occasions : if the thermometer has not been calibrated for several years and shows a drift
Metrological services : verification or calibration in 3 points (+0°C or initial value) or customers's specifications
Method : calibration bench consisting of a thermostatic bath and a traceable chain Pt100 4 wires

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CO₂ DETECTOR

Type : portable CO₂ detector
Recommended frequency : on an annual basis
Metrological services : verification and calibration
Method : calibration in 3 points (ambient air, zero, and air saturated the with reference gas)

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OVENS : 0 TO 1600°C

Type : any kind of oven
Recommended frequency : on an annual basis, or after a maintenance or repair operation
Metrological service : characterization of the oven
Method : norms FDX 15-140 and NF EN 60068-3-11

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HYDROMETERS – ALCOHOLMETERS

Type : glass hydrometers and alcoholmeters
Recommended frequency : according to internal norms of the user institution
Metrological services :

- simplified hydrostatic method, in 3 points, at 15°C or 20°C, uncertainty 0,04 g.dm³ + 0,3 division
- complete hydrostatic method, in 3 points, at 15°C or 20°C, uncertainty 0,02 g.dm³ + 0,1 division if length between 2 marks > 1 mm

Method : Cuckow calibration bench or by standard comparison

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LNE has officially authorized Laboratoires Dujardin-Salleron to carry out initial verifications for hydrometers and alcoholometers.

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PRESSURE : -1 TO +15 BAR

Type : Bourdon tube pressure gauge
Recommended frequency : on an annual basis
Metrological service : calibration in 5 or 10 points, depending on the class and the diameter of the pressure gauge
Method : by standard comparison

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MICROPIPETTES : FROM 200 µL TO 10 ML

Type : piston-operated volumetric apparatus
Recommended frequency : according to internal norms of the user institution
Metrological service : verification in 1 point
Method : gravimetry according to ISO 8655 norm

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
VOLUMETRY : FROM 200 µL TO 10 L

Type : laboratory glassware, graduated or gauged
Recommended frequency : according to internal norms of the user institution
Metrological services :

- control in 1 point with mercury or water
- verification at the mark with mercury or water
- calibration at the mark with mercury or water
- verification in 3 points with mercury or water
- calibration at the mark with mercury or water

Methods :

- graduated cylinders according to NF EN ISO 4788
- 1-mark volumetric flask according to NF EN ISO 1042
- burettes according to NF EN ISO 385
- laboratory glassware according to NF EN ISO 384
- volumetric instruments according to NF EN ISO 4787

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ALPHABETICAL INDEX

A

Acetaldehyde 33
Acetic acid
 Kit for titration (vinegars) 12
 Spectrophotometry 33
ALCO 02 26
Alcohol
 Alcoholmetry 28
 Distillation 23
 Ebulliometry 21
 Pycnometry 30
Alcohol lamp 92
Alcoholmeters
 "chercheur" 28
 class II verified in 3 points 28,29
 in case 29
 precision verified in 3 points 28
 small model 28
 with certificate 29
Alcoholmetric tables in degree 128
Anaerobic culture 38
Anthocyanins 34
Aphrometers 59
Apron 112
Ascorbic acid 33
Aspiration/Oxydation 17
ATP-meter 109
Autoclaves 85

B

Balances 86
Beakers 66
Benches 98
Bentotest 56
Books 129
Botrytest 31
Bottles 117
Box of 12 wicks 21
Brushes 81
Buffer solutions 48
Bulb for auto-zeroing burette 64
Burettes
 auto-zeroing 7,64
 digital 65
 graduated in acetimetric degrees 65
 graduated in mg/l of SO₂ 64
 Mohr 7,64
 Schilling 64

C

Calibration 132
Caps (screw assembly) 80
Carbodoseur 59
Centrifuges 91
Centrifuge tubes 91
Ceramic gauze 92
Checkstab 57,61
Chemical products 123
Chromatography on paper 12
Cider hydrometer 50
Citric acid 33
Cleaning solution 48
Clogging index 57
CO₂ portable detector 110

Conductimeters 46
Conductivity standards 48
Cooling bath with circulation 27
Copper 33
Cover glass 37
Cristallizing dishes 65
Cuves for spectrophotometry 34
Cylinders
 calibrated 67
 graduated 68
 without graduation 67
 with thermometer slot 67

D

DE 2000 25
DE COMPACT Marcs 26
DE EVO 25
Degasing kit 77
Density meters (electronic) 30,52
Dissolved CO₂ 59
Distillers-Extractors
 ALCO 02 26
 DE 2000 25
 DE COMPACT (Marcs) 26
 DE EVO 25
 Mothe style 23
Dispensers 72
Dropping pipette 69
Drying racks 81
DS Œnophoto 33
DS CO₂ 110
DS Iodo 15
DS Sulfi 15
DS Temp 54
Dynamometers 87

E

Earplugs 110
Ebulliolog 22
Ebulliometers
 Traditional ebulliometer 21
 Ebulliometer 21
 Ebulliotronic 22

Electrodes

pH 17, 48
platinum 15,17

Electrolyte solution 48

Elevator 93

Enzymatic and colorimetric kits 35

Erlenmeyer 66

Ethycal 125

Extraction apparatus for V.A.

DE2000 25

DE EVO 25

Duclaux-Gayon style 9

Mothe style 9,23

Vola 2000 10

Extracto-œnometers 51

Eye wash station 112

F

Fehling method 19
Fermentest 19
Ferrometer 31

Filter membranes 77

Filter stand 77

Flask holders 75,92

Flasks

Erlenmeyer 66
fill flask 74
filter flask 77
flat/round bottom 75
volumetric 73

Flow detector 78

Flowmeter for Frantz Paul 17

Folded filters 76

Forceps 74,79,81

Frantz Paul apparatus 17,18

Fume hood 105

Funnels 76

Funnel holder 76

Furniture 98

G

Gas burners 92
Gauging scale 128
Glass micropipette 69
Gloves 112
Gluconic acid 33
Glucose 33
Glycerol 33
Growing media 38
Guyot glucometer 50

H

Heating mantles 93

Hoses 78

Hydrometers

15 grams range 49
30 grams range 49
60 grams range 49
100 grams range 49
Baumé 51
extracto-œnometer 51
for beers 50
for ciders 50
"tirage" 51
with certificate 49

Hygrometers 62

I-J

Incubators 94

Inventory triple booklet 128

Iron

Kit for determination 31
Spectrophotometry 33

L

Labox 131

Laccase activity 31

Lactic acid

Chromatography on Whatman paper 12
Spectrophotometry 33

Laminar flow cabinet 104

LQA Wine 36

ALPHABETICAL INDEX

M

Macropipettes 71
Magnetic stirrers 84
Magnifiers 95
Malic acid
 Chromatography on Whatman paper 12
 Spectrophotometry 33
Masks 110
Master TAL 41
Melli-mustimeter 51
Methanol 62
Methods for ebulliometers 22
Micropipettes 71
Microscopes 95
Mini-colorimeter 31
Mini electric compressor 17
Mini Free Eau EVO 27
Mini Logger USB 62
Mixer (for bags) 85
Mustimeters 50

N-O

Numeration cells 37
Oxymeters 55

P

PAL-wine 41
Parafilm 81
Peracetic acid 108
Peroxides 108
Petri dish 38
pH
 buffer solutions 48
 electrodes 17,48
 Litmus paper 45
 pH meter 45
 pH tester 45
 TITRA EVO 16

Photometer free and total chlorine 108

Photometer total hardness 108

Pipette bulb 70

Pipette controller 70

Pipette rinsing set 70

Pipettes

automatic 71
for musts 69
graduated 69
Pasteur 69
serologic 69
volumetric 69

Pipette stand 70

Pitchers 116

Platinum electrodes 15,17

Practical guide of alcoholmetry 128

Press for samples 115

Pressure 60

Protective faceshield 110

Proteic stability 57

Pump for pipettes 70

Q

Quick filtration 76

R

Rack for alcoholmeters 29
Rack for glasses 81
RCT 0-22 41
Reagents 121
Refractometers
 handheld 41
 for bench 43
Refrigerated baths 90
Residual sugars
 Fermentest 19
 Fehling method 19
 Spectrophotometry 33
Ring for stand 79
Ripper method 17
Rods 79
Ruler for barrel 115

S

Safety cabinets 102
Safety glasses 110
Salinometers 51
Scales 86
Slides 37
SO₂
 Extraction and determination 17,18
 Potentiometric method 15,16
 Ripper method 13
Sorbic acid 12
Spark M20 62
Spatula 81
Spectrophotometers 34
Spectrophotometry 33
Spirit sampling 115
Stainless steel temperature probes 54
Stands 79
Standard samples 125
Standard solutions 120
Standard wine 125
Steel plates 79
Stereomicroscope 95
Stills
 all degrees 23
 plants and liquors 23
 simple 23
Stir bars 84
Stoppers (grey-blue) 78
Stoppers (screw assembly) 80
Storage solution 48

T

Tanks 116
Tartaric acid 34
Tartaric stability 61
Test tubes 74
Test tube rack 74
Thermometers
 digital 54
 enclosed scale 53
 for ebulliometer 21
 infra-red 54
 stem type 53
 thermo-hygrometer 62
 with protective sleeve 53

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for acetic acid in vinegars 12
for Frantz Paul 18
for free and total SO₂ 13
for total acidity 7
for volatile acidity 11

Titravin 124

Total acidity

Kit for titration 7,13
TITRA EVO 8,¹⁶

Trap for CO₂ and acid vapors 10

Tube water-wine for ebulliometer 21

Turbidimeters

portable 56
bench type 57

U-V

Ultrasonic baths 88
Vacuum filtration 77
Vola 2000 10

Volatile acidity

See also Extraction apparatus for V.A.
 Kit for titration 9
 TITRA EVO 8

Vortex stirrer 84

W

Wall cabinets 98
Wash bottles 118
Water baths 89
Water-jet pump 77
Water still 94
Whatman paper 12
Weighing scoops 81
Weight rings 90
Wine thieves 114,115
Wooden tong 74

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

XVIIIème

The first alcoholmeter appeared during the Roman times. Wine, vineyards and measuring instruments have been known for centuries. During the eighteenth and nineteenth centuries, the evolution of knowledge about wine and wine-making took a new dimension.

Chemists (pharmacists or physicians) such as Antoine Baumé (1728-1804) and Joseph, Louis Gay-Lussac (1778-1850) put together the basis of the first biochemical formulas by stating that sugar is converted into alcohol and carbon dioxide.

XIXème

Pasteur (1822-1895) was the first to demonstrate that fermentation is due to microorganisms.

1855

During this same period, Jules Salleron (1829-1897) created a company specialised in precision instrumentation at 1, rue du Pont Lodi in Paris in 1855. He invented his first still, the acetimeter followed by his renowned ebullimeter still used today. In October 1860, the company moved to 24, rue Pavée in the Marais district of Paris, where it remained for more than a century.

24, rue Pavée was the Mansion of the president of the Paris parliament, G. de Lamoignon in 1658. Today, it houses the Historic Library of the City of Paris.

1880

Around 1880, Jules Dujardin (1857-1947) joined the company.

Author of many leaflets and research papers, he took over in 1889 and associated his name to that of Salleron to form the Dujardin-Salleron Establishments. The tandem brought to the Dujardin-Salleron label a scientific aura that their successors have always been very keen on perpetuating.

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1987

The company, now named Laboratoires Dujardin-Salleron, is not a family business any more since 1987, and this is to even better excel in the precision instrumentation field, serving œnology and spirits. The head office moved in Paris (Rue Payenne around 1937, Rue de la Durance from 1987) and in the region of Paris (Arcueil from 1994).

2012

Since 2012, the company is situated in Noizay, in the region of Tours in Loire Valley. This way it is geographically associated to Verrerie Dumas, a glass blowing company, to create synergies around precision instrumentation.

2024

Today, the company has retained its benchmark status in the field of precision instrumentation serving œnology and spirits thanks to maintaining its high standards and its investment in research. Indeed, while continuing its Ecological and Energy Transition, Laboratoires Dujardin-Salleron continues to innovate, in particular thanks to its on-site R&D department ensuring that it anticipates the needs of the sector while guaranteeing you excellent professional equipment.

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Always at your disposal to provide you with the advice you need. Benefit from the experience of our multilingual sales team to guide you in the choice of your equipment.

Service technique / Technical service

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Our team of experienced technicians assembles Dujardin-Salleron analysis devices directly in our workshop to guarantee reliability and precision.

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Keen to support you and offer you the fastest delivery times, we have stock of all the items necessary for your daily analyses.

Service après-vente / After sales service

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We are committed to extending the life of your devices. We carry out maintenance on Dujardin-Salleron equipment on request.

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Our R&D department at the cutting edge of innovation work to offer you professional equipment of excellence and to provide the response to your needs.

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To ensure you control your results at all times, we are supported by the company E.R.IN. for the calibration, control, verification and maintenance of your instruments.

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Our company, which has long been an ambassador for Made in France, has begun its Ecological and Energy Transition in order to decarbonize its activity.

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Laboratoires Dujardin-Salleron are committed in a process of continuous improvement which results in the renewal of ISO 9001 certification.

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Our values are our daily driving force:

- Collaboration: communicate and work as a team towards the same objectives.
- Experience: offer the best service to our customers thanks to our material and living heritage.
- Adaptability: knowing how to learn and evolve according to the professions and the market.

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The content of this catalogue is strictly informative, and is not contractually binding for Laboratoires Dujardin-Salleron. Specifications might change at any time, without notice.

The rights of information, names and registered brands, more generally all documents on this catalogue, are reserved.

The company cannot be held responsible for errors in typography.

Photos and pictures are not contractually binding.

Laboratoires Dujardin-Salleron
872, route de la Gare
37210 NOIZAY - FRANCE
+33 (0)2 47 25 58 25
info@dujardin-salleron.com
www.dujardin-salleron.com

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© Laboratoires Dujardin-Salleron