

DE Compact

Automatic distiller-extractor by steam extraction
Ref. 140850

New generation

Equipped with high dimensions glass parts and a powerful steam generator, the DE Compact carries out the distillation of alcohol and the extraction of volatile acidity according to European regulation and OIV standards, for an unequalled value for money.

Advantages

- **Complying with european regulations** (regulation EEC n°2676/90 and O.J. L272 of 03/10/90) and OIV standards described in the Compendium of International Methods of Analysis of Wines and Musts.
- **Distillation of wines and beverages** which contain until 20% Vol. alcohol
- **Extraction of volatile acids**
- Fast starting and reliable results

Main features

DISTILLATION / EXTRACTION :

- Selection of the time of extraction with a timer
- Automatic stop of the pump supplying water to the steam generator and of the circulating cooling system
- Instantaneous steam production thanks to an integrated steam generator with powerful heating system (2800W)
- High dimension glass parts (rectifying and cooling columns with high efficiency)
- Constant display of the temperature of the refrigeration liquid
- Short work time : 4 to 6 min by operation

SAFETY / MAINTENANCE :

- Protective Plexiglas door in front of the warm glass parts
- Pressure switch to determine low pressure of the refrigeration circuit
- Low maintenance needed, few parts subject to wear and tear
- Easy access to the parts of the device (glass parts, steam generator)

ACCESSORIES INCLUDED :

- Glassware (distillation flask and volumetric flask)
- Connection hoses for the cooling system and for the evacuation of the samples
- Detailed user manual

REQUIRED FACILITIES AND PRODUCTS :

- Power supply : 220V - 50 Hz
- Water supply : minimum pressure 1,2 bar, advisable temperature 18°C, minimum flow 5-6L/mn
- Demineralized water

Principle :

A steam generator fed with demineralized water enables, with steam bubbling into the sample, an extraction of volatile acids, alcohols, etc.

A rectifying column filled with Raschig rings separates the vapours and prevents the distillation of lactic acid. A cooling coil column, equipped with a special tube for connection with the air to equilibrate pressure, condenses and collects all vapours.

The distillate comes out by a glass tip and is collected in a volumetric flask. The extraction stops at the end of a period previously selected on a timer.



Dimensions (HxWxD) and weigh : 95 x 35 x 25 cm ; 20 Kg
Developed and produced in France



LABORATOIRES DUJARDIN-SALLERON



872 route de la Gare
37210 NOIZAY - FRANCE
Tél : +33(0)2 47 25 58 25
Fax : +33(0)2 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com

Entreprise certifiée
ISO 9001 : 2008



FRANCE
Certificat n°2005101402