

Measure of CO₂

Carbodoseur - Van Slyke


The measure of CO₂ can be easily and simply done volumetrically, either with a Carbodoseur, or by mercury suction, according to the Van Slyke method.

Carbodoseur

This method is the simplest and best adapted method for the control of the wine warehouse.

The measuring is done with a plugged cylinder, beforehand filled with the sample to analyze. Once shaken, the cylinder spills out a certain volume of the sample, proportionately with the volume in CO₂.

The volume of wine remaining in the cylinder along with the temperature of the sample, reported on a correlation table, enable to determine the volume in CO₂ with an average precision of 50mg/l.



Ref. 199000

Van Slyke



Ref. 199100

In laboratories, the Van Slyke equipment is used to simply measure the volume in CO₂, even lower than 200mg/l, by mercury suction.

The result is directly read on the burette, graduated in CO₂/liter for a sample of 2 ml.

This method gives very precise measures and a good reproducibility.

The instrument is delivered without mercury.



Mercury is the subject of strict restrictive measures. We cannot be blamed in case of failure to follow the use instructions.

Laboratoires Dujardin-Salleron



872, route de la Gare
37210 NOIZAY - FRANCE
Tel : +33 (0)2 47 25 58 25
export@dujardin-salleron.com
www.dujardin-salleron.com