


Hygrometers



Our hair hygrometers and thermo-hygrometers are accurate instruments to control the temperature and the room humidity. Small-sized, they can be hung, or simply put in a corner of the room you wish to control.




Ref. 100112



Ref. 100110

Specific features

| | Hair hygrometer | ThermoHygrometer | |
|------------------------|-----------------|------------------|--------------|
| T° of use | -35 ... +65°C | -10 ... +60°C | 14 ... 140°F |
| T° scale | | -10 ... +60°C | 14 ... 140°F |
| T° accuracy | +/- 0,5°C | +/- 1°C | +/- 1,8°F |
| Humidity rate scale | 0 ... 100 % | 10 ... 99 % | |
| Humidity rate accuracy | +/- 2 % | +/- 5 % | |
| Dimensions | diameter 10 cm | 110 x 70 x 15 mm | |



LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com