

Refractometer WM7



This digital refractometer measures musts with 7 scales (probable alcohol content, sugar content in Brix, Oeschlé, Baumé, Babo...)

Portable, and very reliable, it has an automatic temperature compensation, a stainless steel sample plate and a liquid-crystal display.



Ref. 196716

Specific features

	WM7		
Scales	Measure range	Resolution	Accuracy
Brix (%)	0 to 45 %	0,1 %	± 0,1 %
T.A. 1990*	0 to 26 %	0,1 %	± 0,2 %
T.A. 1971	0 to 25 %	0,1 %	± 0,2 %
Oe Ger (Allemagne)	0 to 240°	1°	± 1 %
Oe (Oeschle)	0 to 240°	1°	± 1 %
KMW (Babo)	0 to 40°	0,1°	± 0,2 %
Baumé	0 to 21°	0,1°	± 0,2 %
T° measure range	5 to 40°C		
Power supply	1x 9 Volts PP3 battery		
Response time	3 seconds		
Dimensions and weight	190 x 90 x 40 mm - 300g		

* 1% vol. of alcohol corresponds to 16,83 g of sugars per liter of must.

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France

Tél : 02 47 25 58 25

Fax : 02 47 25 58 30

info@dujardin-salleron.com

www.dujardin-salleron.com