

Checker - Checktemp

pH and Temperature tester



Our Check' range

Advantages

- Compact
- Digital display
- Favourably Replace pH litmus paper
- Replacement electrode for the Checker
- Stainless steel probe for checktemp and checktemp 1.



Checktemp 1
Ref. 410102

Checktemp
Ref. 410100

Specific features

Electronic thermometers Checktemp and checktemp 1			
	Measure range	Resolution	Accuracy
T° (C)	-50/+150°C	0,1°	+/-0,3°C
T° (F)	-58/+302°F	0,1°	+/-0,5°F
Power supply	1 button battery 357A (1,5V)		
Dimensions / Weight	66 x 50 x 25 mm / 50g		
T° probe	stainless steel		



Checker
Ref. 701100

Testeur de pH - checker			
	Measure range	Resolution	Accuracy
pH	-50/+150°C	0,1°	+/-0,3°C
Calibration	In 2 points, pH 4.01 / 7.01 or 10.01		
Power supply	2 button batteries 357A (1,5V)		
Dimensions / Weight	66 x 50 x 25 mm / 63g		

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com