

Five Go / Five Easy

pH-meters



Five Go

Ergonomic, with a probe clip to fix to the left or right side of the device, easy-to-use, the pHmeter FiveGo is an equipment adapted to all situations and users.

Equipped with a pH/T° probe, this quality pHmeter meets the accuracy and comfort expectations of professionals.

Delivered with all accessories, ready to use, with probe clip, loop, pH probe with T° sensor, and 2 buffer solution sachets (4.01 and 7.01.)



Ref.708200

Five Easy

Very compact, with a practical probe-holder which can be assembled on the left or right side of the device, easy-to-use, the bench pHmeter FiveEasy is an accurate laboratory device with a good value for money.

This instrument is delivered with all accessories, ready to use, with a probe holder, a network adapter, a pH probe with T° sensor and 2 buffer solution sachets ((4.01 et 7.01).

Advantages

- Data memory for 30 results
- Automatic detection of the buffer during pH calibration
- Automatic or manual temperature compensation
- Simultaneous display of the temperature measure and capture of the final point

Specific features

| | Measure range | Resolution | Accuracy |
|---------------------|---|------------|-----------|
| pH | 0 to 14 | 0,01 | +/- 0,01 |
| mV | -1999 to 1999 | 1 | +/- 1 |
| T° | 0 to 100°C | 0,1°C | +/- 0,5°C |
| Calibration | In 2 points, 3 groups of predefined buffers | | |
| Alimentation | Five Go : 4 AAA batteries Five Easy : 100-240V/50-60 Hz - 9V battery PP3 | | |
| IP Class | IP 54 (waterproof and dust-tight) | | |
| Dimensions / weight | Five Go : 169 x 82 x 36 mm / 180 g Five Easy : 200 x 175 x 52 mm / 600 g | | |

Ref. 758200



LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com