

Macropipettes - Micropipettes

The macropipettes and micropipettes Acura, with digital display, are conceived to meet the needs of the most demanding users.

Ergonomic design for an easy hold and a bigger comfort of use ; guaranteed precision ; solid materials with resistance to shocks, heat, chemicals and UV rays... These pipettes will be quickly indispensable to your laboratory.



Macropipette
Ref. 605052



Micropipette
Ref. 605050

Advantages

- Digital display of the volume
- Pipette filling piston guaranteed without effort
- Adjustable mechanism to eject the tip
- Adjustment of the volume by turning the push button with safety notches to avoid any accidental change in the volume when using the instrument
- Anti-rising protection filters to avoid the pipette contamination
- Instantaneous calibration system under the color cap
- Possibility to put adapters for glass Pasteur pipettes
- Autoclave sterilization without dismantling at 121°C.



Work station
Ref. 605080

Specific features

Ref.	Volume	Division	Tip	Weigh
Macropipettes - Variable volume				
605051	0,5 - 5 ml	0,01 ml	5 ml (605077)	115 g
605052	1 - 10 ml	0,01 ml	10 ml (605078)	130g
Micropipettes - Variable volume				
605053	20 - 200 µl	0,2 µl	200 µl (605075)	90 g
605050	100 - 1000 µl	1 µl	1000 µl (605076)	95 g
Macropipettes - Fixed volume				
605061	2 ml		2 ml	110 g
605062	5 ml		5 ml (605077)	115 g
605063	10 ml		10 ml (605078)	130 g

Optional accessories

- Plastic work station for 3 Acura macro or micropipettes Ref.605080.
- Protection cellulose filter for macropipettes Ref.605070 and 605071.
- Pipette filling tips 200µl, 1000µl, 5ml, 10ml. Ref. above.

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com