

Precision balances



Electronic precision balance with LCD and digital filter.

Specific features

- Stainless steel plate
- Control weight
- Adapter 230 V / 50 Hz
- Protection sleeve on the keyboard and carter
- RS232C output for printer or PC
- Integrated programmes : piece counting, percentage-determination



Characteristics

	601700	601702	601705	601706
Range (g)	200	400	2000	6000
Reading (g)	0,01	0,1	0,1	1
Reproductibility (g)	0,01	0,1	0,1	1
Linearity (g)	± 0,02	± 0,02	± 0,2	± 2
Dimensions plate (mm)	Ø 105	130 x 130	130 x 130	170 x 150
Dimensions (mm)	165 x 230 x 80			
Weight (kg)	1			

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com