

Compact balance



Simplicity, battery life and performances of these compact balances guarantee an excellent value for money.

Original, with exclusive design, and for a universal use, the 2 models of the range will find a place at the production control, in situ control, or anywhere you need a small and portable balance.

Advantages

- Ready-to-use
- Wide display with strong contrast
- Batteries included
- Auto-off function to preserve the batteries
- Tare function facilitates formulation work
- Particular flat design
- Secure and non-slip positioning with rubber feet



Specific features

	601704	601707	601708
Range (g)	2200	500	5200
Reading (g)	1	0,1	1
Repeatability (g)	1	0,1	1
Linearity (g)	± 2	± 0,2	± 3
Available units	G		
Optimal temperature	5-35°C		
Power supply	9 V (included) operating time approx. 12h		
Plastic plate (mm)	diam; 150 mm		
Dimensions (mm)	170 x 240 x 39		
Weight (kg)	0,5 kg		

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com