

Hotplate magnetic stirrer



Advantages

- White Ceramic top with excellent chemical and temperature resistance, easy to clean.
- Cast aluminium body shaped for stability
- «HOT» warning light when the plate temperature is above 50°C (safety for user)
- Spillproof design
- Overheating prevention (independent safety circuit)
- Fast heating and accurate temperature stabilization
- Powerful stirring up to 15 litres
- Fitting for retort rod 12,5 mm (not included)
- Supplied with 2 x 25mm PTFE coated stir bar



Technical features

Ref.100051

	Hotplate magnetic stirring
Plate material	Ceramic
Plate dimensions	160x160 mm
Heated area	120x120 mm
Heater power	500 W
Max. plate temperature	450°C
Stirrer speed	100 to 1500 rpm
Maximum stirring capacity	15 litres
Overall dimensions	190 x 300 x110 mm
Power supply	230 V / 50 Hz / 750 W
Net weight	3,4 Kg

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com