

Seven Go Duo

Professional pH-meter/conductimeter



Easy-to-use, pH/conductivity meter SevenGo Duo is really indispensable for all analysis specialists, and is particularly adapted to food industry. This dual-channel portable instrument measures simultaneously 2 parameters (pH and conductivity).

Advantages

- 3-point pH calibration with 4 pre-defined buffer group, and 1 user-defined buffer group.
- Data storage for 99 GLP data records
- User-friendly instrument thanks to a high-contrast display, with high dimensions digits.
- High level of reproductibility of results thanks to automatic detection of end point, automatic detection of buffers and automatic temperature compensation.
- Steady and waterproof (IP67).

SevenGo Duo pH/Conductivity is delivered in a field kit, including 2 electrodes (pH and conductivity), IP67 protection and 2-meter cable each, calibration sachets for pH and conductivity, sample flasks, a clip for 2 electrodes and the carrying case.



Technical specifications

	Measuring range	Resolution	Accuracy
pH	0.00 to 14.00	0.1	+/- 0,01
mV	-1999.9 à 1999.9	1	+1
Conductivité	0.10 µS/cm ... 500 mS/cm	0.01 ... 1	±0,01%
T°	-5.0 à 105.0°C with ATC	0,1°C	+/- 0,2°C
Calibration	3-point pH calibration		
Power supply	4 AAA batteries 1,5 V or NiMH 1,3 V		
IP Class	IP 67 (water and dustproof)		
Dimensions and weight	220 x 90 x 45 mm / 325 g		

LABORATOIRES DUJARDIN-SALLERON

872 route de la Gare
 37210 NOIZAY - FRANCE
 Tél : 02 47 25 58 25
 Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com