

"Titrette"

Bottle-top digital burette

The bottle-top burette Titrette® enables quick, reliable and highly accurate titrations without the need of big spaces or a power supply at the laboratory, production site or field.

It uses a very innovative system which provides the same accuracy as the class A burettes. Volumes below 20 ml are displayed with 3 decimals.

Advantages

- Large, easy-grip hand wheels, easy to turn for a slow and sensitive drop-wise titration
- Good stability ensured by a compact and light design
- The dispenser automatically detects whether you are filling or titrating by the direction of hand wheel rotation
- Adjustment without tools with Easy Calibration,
- Calibration schedule for the next calibration,
- Save power with the Auto Power Off function,
- Changing between two or three decimal places
- Easy maintenance and cleaning operations

Technical features

«Titrette» 25 ml	25 ml	12.5 ml	2.5 ml
Accuracy ($\pm$ % / $\pm$ μ l)	0.07 / 18	0.14 / 18	0.70 / 18
Coefficient of variation ($\leq$ % / $\leq$ μ l)	0.025 / 6	0.05 / 6	0.25 / 6
Minimum resolution	0.001 ml Above 20 ml: 0.01 ml		
Use conditions	Temperature: +15°C to +40°C (instrument and solutions) Viscosity : up to 500mm ² /s Maximum altitude: 3000m Relative humidity in air: 20 to 90%		
Power supply	2 microbatteries 1.5V		
Dimensions	240x75x75 mm		



Ref. 600131

Delivered with:

- Conformity and quality certificate
- Telescoping filling tube (length 170 - 330 mm)
- Recirculation tube
- 2 microbatteries 1,5 V
- Bottle adapters PP (GL 45/32, GL 45/S 40, GL 32/NS 29/32)
- 2 amber colored light shield inspection windows
- Operating manual.

Options:

- Amber glass flask 2,5 L (ref.600132)
- PC interface for data transmission to the PC

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com