

Beer refractometers

PAL-1 / PAL Plato



PAL-1 / PAL Plato

Digital refractometers

PAL series gives a result within 3 seconds. Amazingly compact, these refractometers fit in a pocket and can be used at any time for outdoor tests. Watertight (IP65), they can be directly rinsed under a water tap. Precise and hardy, they have an automatic temperature compensation, a stainless sample plate and a wide LCD.

PAL-1 : Quick measurement of Brix degree.

PAL Plato : Quick measurement of Plato degree.



PAL-1
Ref. 196711



PAL Plato
Ref. 196709

Specific features

	PAL-1	PAL Plato
Measurement range	0.0 to 53.0 % Brix	0.0 to 30.0 °P
Accuracy	±0.2 %	±0.2 °P
Repeatability	±0.1 %	±0.1 °P
Measurement T°C	10 to 75 °C	
IP Classe	IP 65 (resistant to water and dust)	
Dimensions / weight	5,5 x 3,1 x 10,9 cm - 100g	



LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France
Tél : 02 47 25 58 25
Fax : 02 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com