

# pH testers

## Gondo - PC5



Ref.705200

### Gondo

#### Handheld electronic pH tester

This waterproof pH tester displays both pH and temperature reading simultaneously. It can be automatically calibrated in 2 or 3 points with buffer solution pH 4, 7 and 10. The Gondo is supplied in a carrying case with buffers pH4 and pH7, and 4 batteries 1.5V.

Gondo	Measurement range	Accuracy
pH	-2.00 to 16.00	+/-0.05
T° C	0 to 90.0°C	0.2°C

### PC5

This tester is a practical tool supplied in case with calibration solutions for the analysis, on the field. Backlight, interchangeable combined electrode with high durability, automatic temperature compensation, autoOff... are so advantages that give it a very good value for money.

PC5	Measurement range	Resolution
pH	-2,00 à 16,00	0,01
EC	0 ... 20 mS/cm	auto
T°	0 ... 50°C	0,1



#### LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare  
37210 NOIZAY - France  
Tél : 02 47 25 58 25  
Fax : 02 47 25 58 30  
info@dujardin-salleron.com  
www.dujardin-salleron.com