

Dispensers

Calibrex



The Calibrex™ line includes a large selection of robust dosing instruments with excellent chemical resistance, high performance and simplified maintenance. Intended for the safe and reproducible liquid distribution in volumes ranging from 1 ml to 50 ml.

This dispenser has a PFA coated plunger preventing the crystallization of chemicals. Made of high-tech materials, it provides for broad chemical resistance and is universally adapted for most laboratory reagents.

Advantages

- Friendly volume setting
- Optimized working position : the instrument rotates 360° for adequate positioning on top of the bottle.
- Colour coding identification
- Integrated calibration key
- Long lasting performance stability
- QR coded chemical compatibility
- Quick disassembling, no tool needed.
- Autoclavable at 121°C fully assembled
- Warranty : 2 years

Technical specifications



Ref.	Volume ml	Division ml	Inaccuracy (err. rel. %) (Min.vol / Mid.vol. / Max.vol.)	Imprecision (CV %) (Min.vol / Mid.vol. / Max.vol.)
600161	1 - 5	0,1	< ± 2.0% / < ± 1.3% / < ± 0.6%	< 0.5% / < 0.35% / < 0.1%
600162	1 - 10	0,25	< ± 1.5% / < ± 1.1% / < ± 0.6%	< 0.5% / < 0.35% / < 0.1%
32 mm base thread. Supplied with 300 mm feed tubing, 80 mm delivery jet with stopper, 28, 40, 45 mm bottle neck adapters, QC certificate and operating instructions.				
600163	2,5 - 25	0,5	< ± 1.5% / < ± 1.1% / < ± 0.6%	< 0.5% / < 0.35% / < 0.1%
600164	5 - 50	1,0	< ± 1.5% / < ± 1.1% / < ± 0.6%	< 0.5% / < 0.35% / < 0.1%
45 mm base thread. Supplied with 350 mm feed tubing, 120 mm delivery jet with stopper, 32, 38, 40 mm bottle neck adapters, QC certificate and operating instructions.				

Options and accessories

- Amber or transparent flasks in 250, 500, 1000 and 2500 ml
- Luer membrane filter to avoid contamination of the chemical.
- Extension tubing to dispense within a 60cm distance
- Threaded or groundjoint adapters with different sizes.

LABORATOIRES DUJARDIN-SALLERON

872 Route de la Gare
37210 NOIZAY - France

Tél : 02 47 25 58 25

Fax : 02 47 25 58 30

info@dujardin-salleron.com

www.dujardin-salleron.com