

PC70 Vio

pH-conductivity portable meter

Delivered with a carrying case for easy transport during on-site measurements, the PC70 Vio pH-conductivity meter combines precision and robustness to ensure reliable measurements on-site and in the laboratory.

Advantages

- Ergonomic device, with a protective shock-resistant rubber case.
- Large display.
- Simultaneous display of all measurement information.
- 6 buttons to easily control all the functions.
- Waterproof USB port for data transfer, PC power or 220V power.
- Automatic calibration.
- Automatic temperature compensation.
- Adjustable stability parameter (Low-Med-High) for high accuracy and/or quick measurement.
- Data recording (500 time-stamped measurements).
- Complete kit in a case including a pH electrode and a conductivity probe with integrated thermometer, power supply, PC-Link data acquisition software and USB cable.



Réf. 710300

Technical features

pH-conductivity portable meter pC70 Vio			
	Measurement range	Resolution	Calibration points
pH	- 2,00 à 16,00	0,1,0.01/+0.02	In 1, 2 or 3 points (USA or NIST)
Conductivity	0,00...199,9 mS	0,01/0,1/1 µS 0,01/0,1 mS (automatic)	In 1 and 5 points 84 µS – 1413 µS – 12,88mS – 111,8 mS
T°	0 à 100°C	0,1°C / ±0,5°C	
Temperature compensation	Automatic or manual		
Power supply	3 AA batteries (1,5V) or external power supply (AC / DC adapter supplied) via USB		
IP Class	IP57		
Dimension and weight	Device : 185 × 85 × 45 mm / 400 g		

Laboratoires Dujardin-Salleron



872, route de la Gare
37210 NOIZAY - FRANCE
Tel : +33 (0)2 47 25 58 25
Fax : +33 (0)2 47 25 58 30
export@dujardin-salleron.com
www.dujardin-salleron.com

