

# pH50 Vio

## Bench pHmeter PH50 Vio

Supplied complete for immediate use, the bench pHmeter PH50 Vio is a perfect compromise between design and ease of use.

### Advantages

- Large display with backlighting.
- pH and temperature parameters simultaneously displayed.
- Only 4 keys allow the user to control all functions.
- Automatic calibration in 2 or 3 points.
- Automatic temperature compensation.
- Adjustable stability parameter (Low-Med-High) for high accuracy and/or quick measurement.
- Constant monitoring of the electrode's functional parameters to alert the operator if adjustment is necessary.
- Supplied with pH electrode with integrated temperature sensor and stand for a perfectly stable electrode.



Réf. 758100

### Technical features

pHmeter PH50 Vio		
	Measurement range	Resolution
pH	0,00 à 14,00	0,1 ou 0,01
T°	0 à 100°C	0,1°C / ±0,5°C
Temperature compensation	Automatic or manual	
Calibration	2 or 3 points among series 1,68- 4,01- 6,86- 9,18 ou 1,68- 4,01- 7,00- 10,01	
Power supply	5V adapter 100 ... 240V power supply	
IP Class	IP54	
Dimension and weight	Device : 160 x 140 x 45(h) mm / 380 g	

Laboratoires Dujardin-Salleron

872, route de la Gare  
37210 NOIZAY - FRANCE  
Tel : +33 (0)2 47 25 58 25  
Fax : +33 (0)2 47 25 58 30  
export@dujardin-salleron.com  
www.dujardin-salleron.com



Intertek