

pH 7 Vio

Portable pHmeter

Delivered in carrying case to make easier on-site checkings, the pHmeter «pH7 Vio» combines accuracy and steadiness to ensure reliable measurements outside, at the warehouse or in the laboratory.

Advantages

- Elegant and ergonomic case, with a practical soft rubber protection.
- Large display.
- pH and temperature parameters simultaneously displayed.
- Only 4 keys allow the user to control all functions in an intuitive and easy way.
- Automatic calibration in 2 or 3 points.
- Automatic temperature compensation.
- Adjustable stability parameter (Low-Med-High) for high accuracy and/or quick measurement.
- All the functional parameters of connected electrode, during measurement or calibration, are constantly checked and inform the operator for easy trouble-shooting.
- Complete kit in test including a pH electrode with temperature sensor inside, colored pH buffer solution and accessories.



Réf. 708100

Technical features

Portable pHmeter PH7 Vio		
	Measurement range	Resolution
pH	0.00 to 14.00	0.1 or 0.01
T°	0 to 100°C	0.1°C / ±0.5°C
Automatic temperature compensation	Manual or automatic	
Calibration	In 2 points or 3 points among series : 1.68- 4.01- 6.86- 9.18 or 1.68- 4.01- 7.00- 10.01	
Power supply	3 batteries AA (1.5V)	
IP class	IP57	
Dimensions and weight	Instrument only : 86 × 196 × 33 mm / 295 g	

Laboratoires Dujardin-Salleron



872, route de la Gare
37210 NOIZAY - FRANCE
Tel : +33 (0)2 47 25 58 25
Fax : +33 (0)2 47 25 58 30
info@dujardin-salleron.com
www.dujardin-salleron.com

